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INDICATIONS, DEVICE DESCRIPTION, HANDLING AND STERILITY

2 **STERILE** See product label for method

INDICATIONS FOR USE WITH CEMENT

The LCS® Complete™ – P.F.C.® Sigma™ RP Mobile Bearing Total Knee System is indicated for cemented use in cases of osteoarthritis and rheumatoid arthritis. The RPF and RPS inserts and femoral components are indicated where a higher than normal degree of post-operative flexion is required. The rotating platform prosthesis and modular revision components are indicated for revision of failed knee prostheses.

INDICATIONS FOR USE WITHOUT CEMENT

The porous coated Keeled and Non Keeled M.B.T.™ (Mobile Bearing Tibial) Tray configurations of the LCS Total Knee System are indicated for noncemented use in skeletally mature individuals undergoing primary surgery for reconstructing knees damaged as a result of noninflammatory degenerative joint disease (NIDJD) or either of its composite diagnoses of osteoarthritis and post-traumatic arthritis pathologies. The Rotating Platform device configuration is indicated for use in knees whose anterior and posterior cruciate ligaments are absent or are in such condition as to justify their sacrifice. The P.F.C. Sigma RP Curved bearings when used with the P.F.C. Sigma Cruciate Retaining femoral component can be used in posterior cruciate ligament retaining procedures.

DEVICE DESCRIPTION

The LCS Complete – P.F.C. Sigma RP Mobile Bearing Total Knee System consists of three components intended to replace the femoral, tibial, and patellar articular surfaces of the knee joint. For component usage compatibility see Table 1 and Table 2.

Femoral Components

The LCS Complete femoral components are manufactured from Co-Cr-Mo alloy conforming to ASTM standard F-75. The geometry of the femoral condylar component matches the tibial bearing component in the major load bearing positions of extension and stance phase flexion. The femoral components are available in right and left configurations. The LCS Complete revision femoral components allow for the attachment of femoral augments, sleeves, and ancillary P.F.C. stem extensions to provide additional stability to accommodate for femoral bone loss. The LCS Complete Modular femoral components allow for the attachment of femoral augments and ancillary P.F.C. stem extensions.

Mobile Bearing Tibial Tray (Primary and Revision)

The M.B.T. trays are manufactured from Co-Cr-Mo alloy conforming to ASTM standard F-75. The metal trays have a hollow, conical intramedullary stem. The flat superior surface of the tray has a central hole continuous with the hollow conical medullary stem which receives the conical spike of both the LCS Complete and P.F.C. Sigma RP bearing components. The bearing component is not secured to the metal tray, but rotates axially to maintain congruent contact with the femoral prosthesis up to 40° of flexion. The M.B.T. tray, both primary and revision, is offered in both porous and cemented versions. The M.B.T. tray is offered in Keeled and Non-Keeled versions, while the M.B.T. revision tray is only offered in a Keeled version. The M.B.T. revision tray allows for the attachment of tibial augments, sleeves, and ancillary P.F.C. stem extensions to provide additional stability to accommodate for tibial bone loss. The M.B.T. trays articulate with the LCS Complete, the P.F.C. Sigma RP, and the LPS tibial hinge bearing components. For component usage compatibility see Table 2.

Rotating Platform (RP) Bearings

All RP bearings are manufactured from UHMWPE (the P.F.C. Sigma RPF and LCS Complete RPS inserts are reinforced with ASTM F-90 Co-Cr-Mo alloy pins, and the P.F.C. Sigma RP TC3 and LCS Complete VVC inserts are reinforced with ASTM F-136 Titanium alloy pins). Both the LCS Complete and P.F.C. Sigma RP (including RPF) bearings articulate with the M.B.T. primary and revision trays. For component usage compatibility see Table 1 and Table 2.

Patellar Components

The metal backed LCS mobile bearing patellar components articulate with the LCS femoral components. The LCS and P.F.C. Sigma One Piece Fixed Bearing UHMWPE patellar components articulate with the LCS and P.F.C. Sigma femorals respectively.

FIXATION SURFACE FINISHES

Porous Surface Finish

The entire fixation surfaces of all three metal components are covered with a porous coating of sintered microbeads manufactured from the same Co-Cr-Mo alloy. The porous coating is composed of small, spherically-shaped metal particles metallurgically bonded together and affixed to the prosthesis by a proprietary sintering process. The porous structure is designed to act as a scaffold in which tissue may grow, providing a means for fixation of the device.

Textured Surface Finish for Cemented Application

Cemented LCS Co-Cr-Mo femoral and M.B.T. tray components are available for cemented applications. The cemented components have a grit blasted finish to provide for enhanced cement fixation.

HANDLING

It is important that the size of the implant be properly selected initially, as once an implant has been inserted, cleaning and sterilization techniques cannot ensure removal of all debris from the surfaces, and it should not be reinserted. Proper handling of any implant (porous, cemented or UHMWPE) is mandatory. Any prosthesis should be handled at all times by personnel wearing surgical gloves. Prior to surgical use, a visual inspection of each implant for possible imperfections should be performed routinely. Damage or alterations to the implant may produce stresses and cause defects that could become the focal point for implant failure. Do not allow contact of the prosthesis with hard objects. A prosthesis should never be reimplanted. Even though the implant appears undamaged, it may be fatigued from previous stresses and may have developed microscopic imperfections that could lead to implant failure.

HOW SUPPLIED

All components of the LCS Complete – P.F.C. Sigma RP Mobile Bearing Total Knee System are supplied prepacked and sterile. The integrity of the packaging should be checked to ensure that the sterility of the contents is not compromised. Remove implants from packaging, using aseptic technique, only after the correct size has been determined.

For metal components only: If the sterile implant is determined to be aseptically compromised but still acceptable for intended use based on physician determination, the implant must be rinsed and sterilized prior to implantation according to the following instructions. POLYETHYLENE COMPONENTS SHOULD NOT BE RE-STERILIZED. DO NOT USE POLYETHYLENE COMPONENTS IF THE STERILE PACKAGE APPEARS TO BE DAMAGED.

STERILIZATION (METAL COMPONENTS ONLY)

RINSING/CLEANING

Use sterile room temperature water or physiological saline to soak the implant. Soak the implant for a minimum of 5 minutes. Immediately dry the product. Inspect the implant prior to sterilization.

If sterilization of a metal component is necessary, the following parameters are recommended as they have been validated for a Sterility Assurance Level (SAL) of 10⁻⁶:

Method	Cycle	Temperature	Exposure Time
Steam	Pre-vacuum	270°F(132°C)	10 minutes

NOTE: The adequacy of any sterilization procedure must be suitably tested. It is critical that appropriate process parameters be validated for each facility's sterilization equipment and product/load configuration by persons who have training and expertise in sterilization processes to substantiate the process and its reliability and reproducibility. Flash sterilization may be conducted, if applicable, according to the specific health care facility's policy.

CONTRAINDICATIONS, WARNINGS, PRECAUTIONS AND ADVERSE EFFECTS

CONTRAINDICATIONS FOR USE WITH AND WITHOUT CEMENT

The use of the LCS Complete – P.F.C. Sigma RP Mobile Bearing Total Knee System is contraindicated in:

- the presence of osteomyelitis, pyrogenic infection or other overt infection of the knee joint. Every effort should be made to rule out the possibility of preoperative sepsis in a patient who has one or more of the following abnormalities:
 - fever or local inflammation;
 - rapid destruction or bone resorption apparent on x-rays;
 - elevation of the erythrocyte sedimentation rate or white blood cell count unexplained by other disease or a marked shift in the white blood cell differential count.
- patients with any active infection at sites such as the genitourinary tract, pulmonary system, skin or any other site. Should a patient have any infection prior to implantation, the foci of the infection must be treated prior to, during and after implantation.
- patients with loss of musculature or neuromuscular compromise leading to loss of function in the involved limb or in whom the requirements for its use would affect recommended rehabilitation procedures.
- patients with severe osteoporosis or other metabolic bone diseases of the knee;
- patients with any of the following conditions:
 - lesions of the supporting bone structures (e.g. aneurysmal or simple bone cysts, giant cell tumor or any malignant tumor),
 - systemic and metabolic disorders leading to progressive deterioration of solid bone support,
 - the presence of severe instability secondary to advanced loss of osteochondral structure or the absence of collateral ligament integrity, fixed deformities greater than 60° of flexion, 45° of genu varus or valgus,
 - known drug or alcohol addiction,
 - skeletally immature individuals and the presence of allergic reaction to implant metals or polyethylene are also contraindications for the noncemented, porous coated, M.B.T. and LCS Complete – P.F.C. Sigma RP Mobile Bearing device configurations, and for the cemented use of all device configurations of the LCS Complete – P.F.C. Sigma RP Mobile Bearing Total Knee System.

CONTRAINDICATIONS FOR USE WITHOUT CEMENT

(All Porous Coated Keeled and Non-Keeled M.B.T. Tray Device Configurations Only)

Noncemented use of the Porous Coated Keeled or Non-Keeled M.B.T. Tray device configurations is contraindicated in patients with sufficient loss in quantity or quality of bone stock (as determined on x-ray) such that successful noncemented fixation is unlikely. Additional contraindications may become apparent at the time of surgery.

These include:

- vascular deficiency at the bone site;
- inadequate bone stock to assure both a firm press fit and close apposition of the cut bone surfaces to the prosthesis;
- the inability to make bone cuts so as to assure both correct component position and intimate apposition of bone and prosthetic surfaces;
- inadequate bone quality (e.g. severe osteoporosis) and lack of stability of the implanted components.

In the presence of any of the above conditions, noncemented implantation of the Porous Coated Keeled or Non-Keeled M.B.T. Tray device configurations of the LCS Complete – P.F.C. Sigma RP Mobile Bearing Total Knee System is contraindicated, and the components should be fixed with cement.

WARNINGS

Familiarity with and attention to the surgical technique utilized with this device are imperative for best results. The correct selection as well as the correct seating/placement of the prosthetic implant is extremely important. During surgery, particular attention to tracking of the patella is also required for a successful result. Thus, failure to use the optimum size implant, failure to adequately seat the component adjacent to adequate bone and failure to ensure that the component is stable may result in dislocation, subsidence, fracture or loosening of the components. The proper size selection, the choice of and careful use of the components and the use of trial prostheses are imperative. The P.F.C. stem extensions can only be used with M.B.T. revision trays and LCS Complete Revision and Modular femoral components. LCS Complete – P.F.C. Sigma RP Mobile Bearing Total Knee components, instruments and trial prostheses should not be used together with those of another manufacturer. Because dimensional compatibility cannot be assured, adverse outcomes can result from the use of components from different manufacturers.

The noncemented use of the M.B.T. tray device configurations of the LCS Total Knee and P.F.C. Sigma Total Knee is a technically demanding surgical procedure that requires careful patient selection.

In particular, it is necessary that there be a close bone/prosthesis interface for the components utilized during the operative procedure (see **CONTRAINDICATIONS FOR USE WITHOUT CEMENT** section).

A postoperative management program is vital. It is recommended that the program be modified according to the condition of the patient and the extent of soft tissue and ligament reconstruction. The safety and effectiveness of the cemented use of the LCS Total Knee in patients under 41 years of age have not been established. The safety and effectiveness of the noncemented use of the M.B.T. tray device configurations in patients under 50 years of age have not been established. The safety and effectiveness of the noncemented use of porous coated components has not been established in patients undergoing revision procedures. The safety and effectiveness of the noncemented use of the M.B.T. tray device configurations for indications other than noninflammatory degenerative joint disease (NIDJD) and in bilateral applications have not been established. The implantation of the P.F.C. Sigma RPF insert and femoral component will not in themselves guarantee a high level of post-operative flexion. The degree of post-operative flexion is multi-factorial. These factors include, but are not limited to, surgical technique, patient build, pre-operative flexion and age.

PRECAUTIONS

The surgeon should discuss all physical and psychological limitations inherent to the use of this device with the patient preoperatively. Particular discussion should be directed to the issues of premature weight bearing, activity levels and the necessity for periodic medical follow-up.

Particular attention should be paid to the handling of the components. Contact between the porous coated components and cloth or other fiber releasing materials should be avoided in order to minimize contamination of the porous surfaces with adherent fibers (see **HANDLING** section for further information).

POTENTIAL ADVERSE EFFECTS OF THE DEVICE ON HEALTH

Multicenter clinical trials of the cemented Bicruciate Retaining and Rotating Platform and the uncemented Posterior Cruciate Retaining and Rotating Platform device configurations of the LCS Total Knee showed the complications occurring were of similar incidences as those occurring in comparison groups. The frequencies of the complications reported during the multicenter trial of the cemented Bicruciate Retaining Meniscal Bearing (n=356) and Rotating Platform (n=562) device configurations are as follows:

Bicruciate Retaining Meniscal Bearing LOCAL INTRAOPERATIVE- 3-bearing size (0.8%); 3-bone graft (0.8%); 20-cruciate bridge fracture (5.6%); 3-early infection present at surgery (0.8%); 1-fracture cement (0.3%); 3-implant size (0.8%); 10-lack of bone stock (2.8%); 14-other (3.9%); 4-patella problem (1.1%); 2-rebuilding (0.6%); 8-soft tissue release (2.2%); 7-technical error (2.0%); 6-tibial fracture (1.7%) **LOCAL POSTOPERATIVE-** 2-bearing extrusion (0.6%); 4-femoral loosening (1.1%); 1-femoral/patellar dislocation (0.3%); 1-femoral/patellar subluxation (0.3%); 2-hematoma (0.6%); 13-infection (3.7%); 2-patella fracture (0.6%); 6-phlebitis (1.7%); 1-peroneal palsy (0.3%); 3-pulmonary embolism (0.8%); 11-tibial loosening (3.1%) **Rotating Platform LOCAL INTRAOPERATIVE-** 2-bearing size (0.4%); 5-bone graft (0.9%); 17-cruciate bridge fracture (3.0%); 9-deformity (1.6%); 4-early infection present at surgery (0.7%); 7-fracture cement (1.2%); 4-implant size (0.7%); 27-lack of bone stock (4.8%); 3-laxity (0.5%); 3-loose component (0.5%); 23-other (4.1%); 8-patella

problems (1.4%); 17-rebuilding (3.0%); 30-soft tissue release (5.3%); 3-technical error (0.5%); 2-tibial fracture (0.4%)

LOCAL POSTOPERATIVE- 2-femoral loosening (0.4%); 2-femoral/patellar dislocation (0.4%); 7-femoral/tibial dislocation (1.2%); 5-femoral/tibial subluxation (0.9%); 2-hematoma (0.4%); 23-infection (4.1%); 1-patella fracture (0.2%); 1-patella loosening (0.2%); 4-phlebitis (0.7%); 4-peroneal palsy (0.7%); 3-pulmonary embolism (0.5%); 1-tibial fracture (0.2%); 7-tibial loosening (1.2%). The frequencies of the complications reported during the multicenter trial of the noncemented porous coated **Tricompartmental Rotating Platform** device configuration where n=395 are as follows: **LOCAL INTRAOPERATIVE-** 1-bone fracture (0.4%); 1-division of medial collateral ligament (0.4%); 1-inaccurate bone resection (0.4%); 6-lateral retinacular release (2.4%); 1-medial collateral ligament laxity (0.4%); 1-Pellegrini-Stieda disease (0.4%); 1-poor bone apposition of implant (0.4%); 6-poor bone stock (2.4%) **LOCAL POSTOPERATIVE-** 5-adhesions/periarticular fibrosis (2.0%); 3-bone fracture (1.2%); 2-dislocation/subluxation (0.8%); 2-effusion/hemarthrosis (0.8%); 1-genu varus deformity (0.4%); 4-hematoma (1.6%); 4-infection (1.6%); 1-instability (0.4%); 2-ligament laxity (0.8%); 1-loose component (0.4%); 1-malaligned femoral component (0.4%); 1-pain (0.4%); 1-patella malposition (0.4%); **SYSTEMIC-** 1-bone fracture (0.4%); 4-cancer (1.6%); 2-embolism (0.8%); 1-foot drop (0.4%); 1-hand numbness (0.4%); 1-hiccups (0.4%); 1-ileus (0.4%); 1-joint replacement not in knee (0.4%); 1-myocardial infarction (0.4%); 1-Parkinson's disease (0.4%); 6-phlebitis/thrombophlebitis (2.4%); 1-pneumonia (0.4%); 1-pulmonary problem, unspecified (0.4%); 1-stroke (0.4%); 1-transfusion reaction (0.4%); 1-hematuria (0.4%); 1-urinary retention (0.4%); 1-urinary tract infection (0.4%). The frequencies of the complications reported during the same multicenter trial for the noncemented porous coated **Tricompartmental Posterior Cruciate Retaining device** configuration where n=463 are as follows: **LOCAL INTRAOPERATIVE-** 5-fracture cruciate bridge (1.1%); 1-extreme laxity (0.2%); 1-implant size (0.2%); 1-other (0.2%); 1-patella problems (0.2%); 3-soft tissue release (0.6%); 1-technical difficulty (0.2%) **LOCAL POSTOPERATIVE-** 18-adhesions/periarticular fibrosis (3.9%); 12-bearing problems (extrusion, dislocation, subluxation, failure) (2.6%); 1-broken drain (0.2%); 1-effusion (0.2%); 1-flexion contracture (0.2%); 1-foreign body proximal tibia (0.2%); 1-hematoma (0.2%); 1-incisional neuroma (0.2%); 4-infection (0.9%); 2-ligament laxity (0.4%); 4-loosening (0.9%); 1-muscle weakness (0.2%); 2-nerve palsy (0.4%); 4-pain (0.9%); 1-patella problems (0.2%); 1-saphenous neuroma (0.2%); 1-scar tissue formation (0.2%); 1-tibial/femoral dislocation (0.2%); 6-thrombosis (1.3%); 1-phlebitis (0.2%); 1-patient fall (0.2%) **SYSTEMIC-** 2-anemia (0.4%); 3-cardiovascular accident (0.6%); 1-congestive heart failure (0.2%); 1-cardiovascular problem (0.2%); 1-pulmonary embolism (0.2%); 1-embolism, deep thrombosis lower leg (0.2%); 1-vascular occlusion (0.2%); 1-pneumonia; urinary problem (0.2%); 1-pulmonary problem (0.2%).

The five year patient survival estimates and 95% confidence intervals for the Rotating Platform cemented and cementless tricompartmental, NIDJD population were 91.3% (85.6-97.0) and 84.6% (79.5-89.6) respectively.¹ The five year patient survival estimates and 95% confidence intervals for the Bicruciate Retaining cemented and Posterior Cruciate Retaining cementless tricompartmental, NIDJD population were 86.6% (80.8-92.4) and 88.7% (84.5-93.0) respectively.

The probabilities of the prosthesis not being revised at five years are 88.3% (81.7-95.0) for cemented and 96.6% (94.2-98.9) for cementless Rotating Platform cases and 83.8% (76.4-92.1) and 94.4% (91.1-97.7) for the Bicruciate Retaining cemented and Posterior Cruciate Retaining cementless cases respectively.

INFORMATION

Where further information is desired, please contact your DePuy representative.

CAUTION

Federal Law (USA) restricts this device to sale by or on the order of a physician.

¹ The comparison for the number of deaths between these groups tested as statistically significant (p=.040).

TABLE 1 – Bearing/Femoral Component Compatibility:

		FEMORAL COMPONENT									
		LCS	LCS Complete RPS	LCS Complete Modular	LCS Complete Revision	LCS Complete	P.F.C. Sigma C/R	P.F.C. Sigma C/S	P.F.C. Sigma TC3	P.F.C. Sigma RPF	S-ROM and LPS
BEARING	LCS Complete VVC				X						
	LCS Complete RP	X	X	X	X	X					
	LCS Complete RPS		X		X						
	LCS Meniscal Bearings	X		X	X	X					
	LCS C/S All-Poly Tibia	X	X	X	X	X					
	P.F.C. Sigma RP Curved						X				
	P.F.C. Sigma RP Stabilized							X	X		
	P.F.C. Sigma RP TC3								X		
	P.F.C. Sigma RPF									X	
	LPS Tibial Hinge Bearing										X

TABLE 2 – Bearing/Tibial Component Compatibility:

		TIBIAL COMPONENT	
		M.B.T. Tray	M.B.T. Revision Tray
BEARING	LCS Complete VVC		X
	LCS Complete RP	X	X
	LCS Complete RPS	X	X
	P.F.C. Sigma RP Curved	X	X
	P.F.C. Sigma RP Stabilized	X	X
	P.F.C. Sigma RP TC3		X
	P.F.C. Sigma RPF	X	X
	LPS Tibial Hinge Bearing		X

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DePuy Orthopaedics, Inc.
DePuy International Ltd