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2 **STERILE** See product label for method**DESCRIPTION**

The DePuy LPS™ Limb Preservation System is designed for the replacement of the mid-shaft portion of the femur, proximal, distal and/or total femur, and proximal tibia. The DePuy LPS system offers a variety of component options (including, but not limited to, proximal femoral bodies, segmental components, distal femoral components, femoral stems, tibial stems, proximal tibial components, hinged tibial insert bearings, metaphyseal sleeves, diaphyseal sleeves, and adapters). The components, which can be used in conjunction with certain components from other systems, are for treatment of patients presenting bone loss and deformity associated with bone tumors resection, trauma, infection, and difficult revision arthroplasty. A total femoral replacement is possible in those cases where no part of the femur can be salvaged.

INDICATIONS

The DePuy LPS is intended for use in replacement of the mid-shaft portion of the femur, proximal, distal and/or total femur, and proximal tibia, especially in cases that require extensive resection and replacement. Specific diagnostic indications for use include:

- malignant tumors (e.g., osteosarcomas, chondrosarcomas, giant cell tumors, bone tumors) requiring extensive resection and replacement;
- patient conditions of noninflammatory degenerative joint disease (NIDJD), e.g. avascular necrosis, osteoarthritis, and inflammatory joint disease (IJD), e.g., rheumatoid arthritis, requiring extensive resection and replacement;
- revision cases for a failed previous prosthesis requiring extensive resection and replacement;
- severe trauma requiring extensive resection and replacement.

The DePuy LPS is also intended for use in bone loss post-infection, where the surgeon has elected to excise the bone and replacement is required.

The S-ROM tibial tray and the non-porous coated straight and bowed stems are intended for cemented use only.

The porous-coated metaphyseal and diaphyseal sleeves are intended for either cemented or cementless applications.

CONTRAINDICATIONS

1. Active local or systemic infection.
2. Loss of bone or musculature, osteoporosis, neuromuscular compromise or vascular deficiency in the affected limb in sufficient degree to render the procedure unjustifiable (e.g., absence of musculoligamentous supporting structures, joint neuropathy).

NOTE: Diabetes, at present, has not been established as a contraindication. However, because of the increased risk for complications such as infection, slow wound-healing, etc., the physician should carefully consider the advisability of knee replacement in the severely diabetic patient.

WARNINGS AND PRECAUTIONS

CAUTION: Federal Law (USA) restricts this device to sale by or on the order of a physician.

- Implants and trial components from different manufacturers or implant systems should never be used together.
- Prosthesis components should never be reimplanted. Even though the implant appears undamaged, the implant may have developed microscopic imperfections, which could lead to failure.
- Always use a trial prosthesis for trial purposes. Trials should not be assembled with any components intended for permanent implantation. Trials must have the same configuration size, etc. as the corresponding components to be permanently implanted.
- Do not alter or modify implants in any way.

CAUTION: The following conditions, singularly or concurrently, tend to impose severe loading on the affected extremity thereby placing the patient at higher risk of failure of the prosthesis.

1. Obesity or excessive patient weight.
2. Manual Labor.
3. Active sports participation.
4. High levels of patient activity.
5. Likelihood of falls.
6. Alcohol or drug addiction.
7. Other disabilities, as appropriate.

The following physical conditions, singularly or concurrently, tend to adversely affect the fixation of prosthesis:

1. Marked osteoporosis or poor bone stock.
2. Metabolic disorders or systemic pharmacological treatments leading to progressive deterioration of solid bone support for the implant (e.g., diabetes mellitus, steroid therapies, immunosuppressive therapies, etc.).
3. History of general or local infections.
4. Severe deformities leading to impaired fixation or improper positioning of the implant.
5. Tumors of the supporting bone structures.
6. Allergic reactions to implant materials (e.g., bone cement, metal, polyethylene).
7. Tissue reactions to implant corrosion or implant wear debris.
8. Disabilities of other joints (i.e., hips or ankles)

A higher incidence of implant failure has been reported in paraplegics and in patients with cerebral palsy or Parkinson Disease. WHEN THE SURGEON DETERMINES THAT A PROSTHESIS IS THE BEST MEDICAL OPTION AVAILABLE AND DECIDES TO USE THIS PROSTHESIS IN A PATIENT WHO HAS ANY OF THE ABOVE CONDITIONS OR WHO IS SIMPLY YOUNG AND ACTIVE, IT IS IMPERATIVE THAT THE PATIENT BE INSTRUCTED ABOUT THE STRENGTH LIMITATIONS OF THE MATERIALS USED IN THE DEVICE AND FOR FIXATION AND THE RESULTANT NEED TO SUBSTANTIALLY REDUCE OR ELIMINATE ANY OF THE ABOVE CONDITIONS.

The surgical and postoperative management of the patient must be carried out with due consideration for all existing conditions. Mental attitudes or disorders resulting in a

patient's failure to adhere to the surgeon's orders may delay postoperative recovery and/or increase the risk of adverse effects including implant or implant fixation failure. Excessive physical activity or trauma to the replaced joint may contribute to premature failure of the prostheses by causing a change in position, fracture, and/or wear of the implants. The functional life expectancy of prosthetic implants is, at present, not clearly established. The patient should be informed that factors such as weight and activity levels may significantly affect wear.

INFORMATION FOR USE

Preoperative

THE SURGEON SHOULD DISCUSS ALL PHYSICAL AND MENTAL LIMITATIONS PARTICULAR TO THE PATIENT AND ALL ASPECTS OF THE SURGERY AND THE PROSTHESES WITH THE PATIENT BEFORE SURGERY.

The discussion should include the limitations and possible consequences of joint replacement, and the necessity to follow the surgeon's instructions postoperatively, particularly in regard to patient activity and weight.

The preoperative planning and surgical techniques for implantation of these prosthetic components evolved from the surgical experience gained during the development of many prostheses. Surgeons should not begin the clinical use of any prosthesis before they have thoroughly familiarized themselves with its specific implantation technique. Certain methods may change with time as further clinical experience is gained. Critical appraisals of such changes are presented at regularly scheduled surgical instruction courses for which periodic attendance is advised. Surgical technique brochures and videos are available from DePuy.

Intraoperative

It is recommended that extra implants be available at time of surgery, including sizes larger and smaller than those expected to be used.

Proper handling of implants is mandatory. These prostheses should only be handled by personnel wearing sterile surgical gloves. Do not allow contact of any component with hard objects. Do not allow porous surfaces to come in contact with cloth or other fiber releasing materials. Prior to use, a visual inspection of each implant component for possible imperfections should be routinely performed. Damage or alterations to any implant component may produce stresses and/or cause defects, which could become the focal point for implant failure.

The use of certain special surgical instruments is suggested in the performance of this surgery. Review of the use and handling of these instruments is important for obtaining prosthesis fixation. Bone excision should be limited to the amount necessary to accommodate the implants. Excessive bone removal may result in mechanical disturbances and bone resorption with subsequent failure of the procedure due to loosening or deformation of the implant. When preparing bone surfaces and positioning the components, proper alignment must be ensured. Patella must be resurfaced when performing a distal femoral replacement.

Prior to closure, the surgical site should be thoroughly cleansed of bone chips, ectopic bone, bone cement, etc. Foreign particles at the metal/plastic articular interface may cause excessive wear. Ectopic bone and/or bone spurs may lead to dislocation or painful and restricted motion. Range of motion should be thoroughly

checked for improper mating, instability or impingement and corrected as appropriate.

Postoperative

Strict adherence by the patient to the surgeon's instructions and warnings is extremely important. Accepted practices should be followed in postoperative care. Postoperatively, the patient should be cautioned to govern his/her activity level to protect the replaced joint from unreasonable stresses. The patient should be released from the hospital with complete written instructions and warnings regarding exercises and therapies and any limitations on his/her activities. Periodic follow-up including x-rays are recommended for close comparison with previous post-op conditions to detect long-term evidence of changes in position, loosening, bending, or cracking of components. If evidence of one or more of these conditions is discovered, patients should be closely observed, the possibilities of further deterioration evaluated, and the benefits of early revision considered.

ADVERSE EVENTS AND COMPLICATIONS

The following are generally the most frequently encountered adverse events and complications in prosthetic implants.

General

1. Early or late loosening, bending, cracking, fracture, deformation, or wear of one or more of the prosthetic components, often related to factors listed under warnings and precautions. Loosening may also occur due to improper fixation or positioning.
2. Fatigue fracture of the implant, often related to factors listed in WARNINGS AND PRECAUTIONS.
3. Early or late infections which may require removal of the implant and a subsequent arthrodesis.
4. Pain, dislocation, subluxation, flexion contracture, decreased range of motion, or lengthening or shortening of the leg caused by improper positioning, looseness or wear of components.
5. Excessive wear of the polyethylene components due to intraoperative damage to the femoral component, loose cement and/or bone fragments and/or high patient activity levels or weight.
6. Fractures of the tibia or femur. Intraoperative fractures are usually associated with revision surgery, deformity and/or severe osteoporosis. Postoperative fractures are usually stress fractures. Fractures can be the result of defects in the cortex due to prior screw holes, misdirected reaming, and/or inadequate or maldistribution of bone cement.
7. Cardiovascular disorders and thromboembolic disease, including venous thrombosis, and pulmonary embolus, and heart attack.
8. Tissue reactions, osteolysis, and/or implant loosening caused by metallic corrosion, allergy or wear debris or loose cement particles.
9. Myositis ossificans, especially in males with hypertrophic arthritis, limited preoperative range of motion and/or previous myositis. The incidence of myositis ossificans is increased with a history of prior surgery and in cases of infection.
10. Possible loss of limb if complications occur.

Early Postoperative

1. Hematoma.
2. Delayed wound healing or wound dehiscence.
3. Varus-valgus deformity.

Late Postoperative

1. Trochanteric avulsion from excessive muscular tension, weight-bearing, or inadvertent intraoperative weakening of the trochanter.
2. Aggravation of problems in the ipsilateral or contralateral knee and ankle joints due to leg length discrepancy, femoral medialization and/or muscular deficiencies.
3. Femoral or acetabular fracture due to trauma or excessive loading, particularly in the presence of poor bone stock caused by severe osteoporosis, bone defects from previous surgery, intraoperative reaming procedures, or bone resorption.
4. Bone resorption, which may contribute to the deterioration of fixation and eventual loosening of the implant.
5. Periarticular calcification or ossification which may lead to a decrease in joint mobility and range of motion.
6. Traumatic arthrosis of the ipsilateral knee secondary to intraoperative positioning of the extremity during surgery.
7. Subluxation or dislocation.
8. Inadequate range of motion due to improper selection or positioning of components, impingement and/or periarticular calcification.
9. Periarticular calcification or ossification, with or without impediment to joint mobility.
10. Patellar fracture as a result of excess tension or inadvertent intraoperative weakening.
11. Aggravated problems of the affected limb or contralateral extremity caused by leg length discrepancy.

The incidence and severity of complications in prosthetic implants are usually greater in revisions than primary operations. Common problems include placement of incision and lack of bone stock. Increased operative time and increased incidence of infection, pulmonary embolism and wound hematoma can be expected with revision procedures.

HOW SUPPLIED

The prosthetic components are individually packaged and supplied STERILE. Remove from the package using accepted aseptic technique only after the correct size has been determined. **Polyethylene components should not be re-sterilized. For polyethylene components: DO NOT USE IF THE STERILE PACKAGE APPEARS TO BE DAMAGED.**

STERILE: Remove from the package using accepted sterile technique only after the correct size has been determined.

CLEANING

For metal components only: If the sterile implant is determined to be aseptically compromised or sterile packaging appears damaged but still acceptable for intended use based on physician determination, the implant must be rinsed and sterilized prior to implantation according to the following instructions.

Use sterile room temperature water or physiological saline to soak the implant. Soak the implant for a minimum of 5 minutes. Immediately dry the product. Inspect the implant prior to sterilization.

STERILIZATION (metal components only)

If sterilization of a metal component is necessary, the following parameters are recommended as they have been validated for a Sterility Assurance Level (SAL) of 10⁻⁶:

Method	Cycle	Temperature	Exposure Time
Steam	Pre-vacuum	270° (132°C)	10 minutes

NOTE: The adequacy of any sterilization procedure must be suitably tested. It is critical that appropriate process parameters be validated for each facility's sterilization equipment and product/load configuration by persons who have training and expertise in sterilization processes to substantiate the process and its reliability and reproducibility. Flash sterilization may be conducted, if applicable, according to the specific health care facility's policy.

CAUTION: Do not re-sterilize polyethylene components.

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DePuy Orthopaedics, Inc.
DePuy International Ltd