IMPORTANT INFORMATION – Please Read Before Use
INFORMATIONS IMPORTANTES - À lire attentivement avant d'utiliser
WICHTIGE INFORMATIONEN - Vor Gebrauch bitte lesen
INFORMACIÓN IMPORTANTE - Sírvase leerla antes de emplear
INFORMAZIONI IMPORTANTI - Leggere prima dell'uso
INFORMAÇÃO IMPORTANTE - É favor ler antes de utilizar
BELANGRIJKE INFORMATIE - Vóór gebruik lezen
VIGTIG INFORMATION - Læses før brug
VIKTIGT! - Läs detta innan produkten används
TÄRKEÄÄ - Lue ennen käyttöä
ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ - Παρακαλούμε διαβάστε πριν από τη χρήση
DŮLEŽITÉ INFORMACE - Před použitím si prosím přečtete tento leták
FONTOS INFORMÁCIÓ - Használat előtt elolvasandó
WAŻNE INFORMACJE - Prosimy przeczytać przed użyciem
DÔLEŽITÉ INFORMÁCIE - Prečítajte si láskavo ešte pred použitím

P.F.C.® Σ and P.F.C. Modular Knee System—Patella
Système de genou P.F.C.® Σ et P.F.C. modulaire—Rotule
P.F.C.® Σ und P.F.C. Modular-Knie-System – Patella
Sistema Modular de Rodilla P.F.C.® Σ y P.F.C. (Rótula)
Sistema protesico di ginocchio modulare P.F.C.® Σ e P.F.C.—Rotula
Prótese Modular do Joelho P.F.C.® Σ e P.F.C. - Rótula
P.F.C.® Σ en P.F.C. Modulair kniesysteem—Patella
P.F.C.® Σ og P.F.C. Modulært knæsystem—Patella
P.F.C.® Σ och P.F.C. modulärt knäsystem—Patella
P.F.C.® Σ ja P.F.C. modulaarinen polvijärjestelmä – Polvilumpio
Σύνθετο Σύστημα Γόνατος—Επιγονατίδας P.F.C.® Σ και P.F.C.
Patelární komponenty P.F.C.® Σ a P.F.C. modulárního kolenního systému
A P.F.C.® Σ és P.F.C. moduláris térd rendszer – patella
System modularny plastyki stawu kolanowego P.F.C.® Σ i P.F.C. – Protezy
Patelárne komponenty – modulárny kolenný systém P.F.C.® Σ a P.F.C.

Contents		
English	Page	4
Contenus		
Français	Page	6
Inhalt		
Deutsch	Seite	8
Contenido		
Español	Página	10
Indice		
Italiano	Pagina	12
Índice		
Português	Página	14
Inhoud		
Nederlands	Pagina	16
Inhold		
Dansk	Side	18
Innehåll		
Svenska	Sida	20
Sisältö		
Suomi	Sivu	22
Περιεχόμενα		
Ελληνικά	Σελίδα	24
Obsah		
Česky	Str.	26
Tartalomjegyzék		
Magyar	Oldal	28
Spis treści		
Polski	Str.	30
Obsah		
Slovenčina	Str.	32
Interpretation of symbols	Page	34
Signification des symboles	Page	34
Bedeutung der Symbole	Seite	34
Interpretación de símbolos	Página	34
Interpretazione dei simboli	Pagina	34
Interpretação de símbolos	Página	34
Betekenis van de symbolen	Pagina	34
Betydning af symboler	Side	34
Förklaring av olika symbolers betydelse	Sida	34
Symbolien selitykset	Sivu	34
Ερμηνεία συμβόλων	Σελίδα	34
Výsvetlivky výrazů a symbolů	Str.	34
Magyarázat a használt jelekhez és kifejezésekhez	Oldal	34
Znaczenie symboli i terminów	Str.	34
Vysvetlenie termínov a symbolov	Str.	34

**STERILE**See product
label for method**Description**

The P.F.C. Σ and P.F.C. Modular Knee System — Patella components are ultrahigh molecular weight polyethylene (UHMWPE) patellar resurfacing prostheses intended for use in conjunction with femoral knee prosthesis components distributed by representatives of DePuy. The articular surface is shaped to match the articular surface of the patellar groove of the femoral components. The P.F.C. Σ and P.F.C. Modular Knee System— Patella components are designed for use with polymethylmethacrylate (PMMA) bone cement.

Indications and Usage

The P.F.C. Σ and P.F.C. Modular Knee System— Patella components are indicated for use in knee replacements for patients suffering from severe pain and disability due to permanent structural damage in the knee joint from rheumatoid arthritis, osteoarthritis, post-traumatic arthritis, collagen disorders, or pseudogout. This damage may also be the result of trauma or failed prior surgical intervention.

The P.F.C. Σ and P.F.C. Modular Knee System— Patella components are indicated for use only with PMMA bone cement.

Information for Use

An instrumentation system, as well as a system of trial components, is available to assure proper fit and alignment of the prostheses. Correct fit and alignment will reduce stresses at interface surfaces to enhance implant fixation. The surgeon should refer to the appropriate surgical technique manual for details on the use of the instrument system and implantation of the prostheses. This manual is available from your local DePuy sales representative or distributor.

Contraindications

Use is contraindicated in cases with active or recent joint sepsis, insufficient bone stock, skeletal immaturity, inadequate soft tissue coverage or support, Charcot arthropathy, significant muscular or neuromuscular compromise or psychosocial disorders that would limit rehabilitation.

Warnings

Improper prosthesis selection or alignment, inadequate fixation, use where contraindicated or in patients where medical, physical, mental, or occupational conditions will likely result in extreme stresses to the implant, may result in premature failure due to loosening, fracture, or wear. Postoperative care is extremely important. The patient should be instructed on the limitations of the device and should be cautioned regarding load bearing, ranges of motion, and activity levels permissible. Early motion and load bearing should be carefully controlled.

Use of other manufacturers' components with this implant is not advised. Use of components other than those recommended could lead to loosening, wear, fracture, and premature failure.

CAUTION: FEDERAL (USA) LAW RESTRICTS THIS DEVICE TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

Precautions

An implant should never be reused. Any implant, once used, should be discarded. Even though it appears undamaged, it may have small defects and internal stress patterns that may lead to failure. Likewise, a new implant should be handled carefully to avoid damage that could compromise the mechanical integrity of the device and cause early failure or loosening.

The wear rate of prosthesis contact surfaces is greatly accelerated if loose fragments of bone cement become detached and act as an abrasive in the bearing surfaces. When using bone cement, care should be taken to remove all excess cement from the periphery of the implant.

The P.F.C. Dome Patella is designed for use with the DePuy femoral components which accept a 1-inch radius of curvature round dome patellar prosthesis.

Do not use the P.F.C. dome patella with the P.F.C. modular knee system cruciate retaining femoral component.

The P.F.C. oval patella may only be used with the P.F.C. Cruciate Retaining femoral component.

The 32mm (extra small) P.F.C. oval patella may only be used with size 1 P.F.C. Cruciate Retaining femoral components.

Adverse Effects

Fracture may occur due to improper preparation of the implant site or if excessive force is used during seating of the implant.

Transient peroneal palsy has been reported following total knee arthroplasty especially after correction of severe flexion or valgus deformities.

Patients have complained of persistent pain and stiffness following total knee arthroplasty. In addition, patellar tendon rupture, femoral-tibial subluxation or dislocation, and persistent ligamentary laxity have been reported with the use of total knee implants. Infection and loosening have been reported following total joint arthroplasty, as have wear and failure due to fracture of knee prosthesis components. Histological reactions have been reported as an apparent response to exposure to a foreign material. The actual clinical significance of these reactions is unknown.

Serious adverse effects may necessitate surgical intervention.

Sterility and Handling**Do Not Reuse**

The P.F.C. Σ and P.F.C. Modular Knee System — Patella components are individually packaged and supplied STERILE. Remove from the package using accepted aseptic technique only after the correct size has been determined.

DO NOT RESTERILIZE. Handle Carefully—Protect from Damage and Contamination.

UHMWPE components must not be resterilized by the hospital because of the possibility of damaging or deforming the plastic.

The operating room staff should be conservative in component selection since opened packages will not be accepted for credit.

How Supplied

The P.F.C. Σ and P.F.C. Modular Knee System — Patella components are supplied in a range of sizes. Please refer to the current price list for product codes and sizes available.

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DePuy Orthopaedics, Inc.
DePuy International Ltd.