IMPORTANT INFORMATION – Please Read Before Use
INFORMATIONS IMPORTANTES - À lire attentivement avant d'utiliser
WICHTIGE INFORMATIONEN - Vor Gebrauch bitte lesen
INFORMACIÓN IMPORTANTE - Sírvase leerla antes de emplear
INFORMAZIONI IMPORTANTI - Leggere prima dell'uso
INFORMAÇÃO IMPORTANTE - É favor ler antes de utilizar
BELANGRIJKE INFORMATIE - Vóór gebruik lezen
VIGTIG INFORMATION - Læses før brug
VIKTIGT! - Läs detta innan produkten används
TÄRKEÄÄ - Lue ennen käyttöä
ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ - Παρακαλούμε διαβάστε πριν από τη χρήση
DŮLEŽITÉ INFORMACE - Před použitím si prosím přečtete tento leták
FONTOS INFORMÁCIÓ - Használat előtt elolvasandó
WAŻNE INFORMACJE - Prosimy przeczytać przed użyciem
DÔLEŽITÉ INFORMÁCIE - Prečítajte si láskavo ešte pred použitím

Sterile **R**

 **DePuy**
a *Johnson & Johnson* company

0086

IFU-0902-00-787
Rev. C, 12/04

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**Description**

The NOILES Rotating Hinge (NRH™) Knee System is a tricompartmental total knee replacement for both primary and revision cases. The NOILES Rotating Hinge Knee System includes the femoral component with hinge pin, the tibial plateau assembly, and the distal femoral augmentation blocks. Replacement bumpers for the femoral component assembly and replacement hinge bearings for the tibial plateau assembly are also available.

The NRH Knee components are designed to be used only with the following S-ROM Modular Total Knee components: femoral sleeve, femoral stem extension, patellar dome, tibial base, tibial sleeve, tibial stem extension, and tibial augmentation blocks. (The femoral sleeve, femoral stem extension, and patellar dome may also be used with the NOILES Posterior Stabilized Knee System.)

The NOILES Rotating Hinge Knee and S-ROM Modular Total Knee Components are indicated for use with polymethylmethacrylate (PMMA) bone cement.

The metal portions of the femoral component and the tibial bases are made of cobalt-chromium (Co-Cr) alloy or titanium alloy (sizes 14 or 24 mm only). The distal femoral augmentation blocks, tibial augmentation blocks, femoral sleeves, tibial sleeves, femoral stem extensions, and tibial stem extensions are made of titanium alloy. The metal portion of the tibial plateau assembly is available in either cobalt-chromium alloy or titanium alloy. All sizes of the tibial sleeves are porous coated; all sizes of the femoral sleeve, except the 20 mm size, are porous coated. The material used for the patellar dome, the bumpers of the femoral component, the plastic portions of the tibial plateau, and the hinge bearings is Ultra High Molecular Weight Polyethylene (UHMWPE). In addition, the patellar dome has an x-ray detectable wire of cobalt-chromium alloy.

Indications and Usage

The NOILES Rotating Hinge Knee is indicated for use with PMMA bone cement in primary or revision cases in patients:

- who have reached skeletal maturity and
- for whom the surgeon has decided to resect both cruciate ligaments or whose cruciate ligaments are absent or incompetent and
- who exhibit insufficiency of lateral/collateral ligaments and other soft supporting tissue due to the following conditions:
 - Rheumatoid arthritis, osteoarthritis, post-traumatic arthritis, and arthritis secondary to a variety of diseases and anomalies
 - Failure of a previous knee reconstruction procedure
 - Trauma

Information for Use

The surgeon should refer to the surgical technique manual for details of assembly technique. An instrument system, as well as a system of trial components, is available to assure proper fit and alignment of the prosthesis. Correct fit and alignment will reduce stresses at interface surfaces to enhance implant fixation. The surgeon should refer to the appropriate surgical technique manual for details on the use of the instrument system and implantation of the prosthesis. This manual is available from your local DePuy Orthopaedics, Inc. sales

representative or distributor.

For porous-coated components, it is recommended that the cement covering the porous coating be used in the low viscosity (i.e., predough) state in order to maximize cement intrusion into the porous surface.

Optimal fixation and implant stability are achieved by maximizing bone coverage. Components are provided in a variety of sizes. The largest available components should be chosen that cover, but do not overhang the femur, tibia, or patella.

Contraindications

1. Active infection or history of general infections or local infectious disease.
2. Vascular insufficiency, muscular atrophy or neuromuscular disease in the affected limb.
3. Advanced loss of osteochondral structure that would preclude proper fixation of the prosthesis.
4. Tumors of the supporting bone structure, systemic and metabolic disorders leading to progressive deterioration of solid bone support.
5. Drug or alcohol addiction, or limiting neuropathic disease.
6. Skeletal immaturity.
7. Obesity or very active lifestyle that can produce loads on the prosthesis that can lead to failure of the fixation of the device or device itself.
8. Allergic reaction to the implant materials.
9. Inadequate flexor and extensor mechanism necessary to achieve a functional prosthetic joint.

Warnings

Improper prosthesis selection or alignment, inadequate fixation, use where contraindicated or in patients whose medical, physical, mental, or occupational conditions will likely result in extreme stresses to the implant, may result in premature failure due to loosening, fracture, or wear. Postoperative care is extremely important. The patient should be instructed on the limitations of the device and should be cautioned regarding load bearing, ranges of motion, and activity levels permissible. Early motion and load bearing should be carefully controlled. The S-ROM tibial base, tibial sleeve, tibial stem extension, and tibial augmentation blocks may not be used with the NOILES Posterior Stabilized Knee.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician.

Precautions

The NOILES Rotating Hinge Knee is designed to articulate from 6° hyperextension to 110° flexion. If, due to grossly inadequate soft tissue integrity, flexion beyond 90° causes luxation of the plateau assembly out of the tibial base, the patient must have a knee brace postoperatively to limit flexion to 90°. In such cases, the surgeon should consider closing the wound with the knee in full extension.

The size of the tibial plateau assembly must correspond with the size of the femoral component.

The size of the tibial augmentation block must correspond to the size of the tibial base.

Femoral sleeves are required when using femoral stem extensions.

A femoral plug is required with the femoral sleeve when a femoral stem extension is not used.

A tibial cap is required with the tibial sleeve when a tibial stem extension is not used.

Tibial augmentation blocks cannot be used when tibial sleeves are being used.

An implant should never be reused. Any implant, once used, should be discarded. Even though it appears undamaged, it may have small defects and/or internal stress patterns that

may lead to failure. Likewise, a new implant should be handled carefully to avoid damage that could compromise the integrity of the device and cause early failure or loosening. The wear rate of prosthesis contact surfaces is greatly accelerated if loose fragments of bone cement become detached and act as an abrasive in the bearing surfaces. When using bone cement, care should be taken to remove all excess cement from the periphery of the implant.

Adverse Effects

Fracture may occur due to improper preparation of the implant site or if excessive force is used during seating of the implant.

Transient peroneal palsy has been reported following total knee arthroplasty, especially after correction of severe flexion or valgus deformities.

Patients have complained of persistent pain and stiffness following total knee arthroplasty. In addition, patellar tendon rupture, femoral-tibial subluxation or dislocation, and persistent ligamentary laxity have been reported with the use of total knee implants. Infection and loosening have been reported following total joint arthroplasty, as have wear and failure due to fracture of knee prosthesis components.

Histological reactions have been reported as an apparent response to exposure to a foreign material. The actual clinical significance of these reactions is unknown.

Serious adverse side effects may necessitate surgical intervention.

Sterility and Handling

Do Not Reuse

The NOILES Rotating Hinge Knee components are individually packaged and supplied **STERILE**. All metal components are sterilized using radiation. Polyethylene components may be sterilized with gas plasma or radiation, as indicated on the outer package label. Remove from the package using accepted aseptic technique only after the correct size has been determined.

Do not resterilize products that are HA-coated or porous-coated; products that are ceramic; or, products that contain plastic components. Resterilization can cause changes to the mechanical and physical properties of these components.

Polyethylene components should not be re-sterilized. For Polyethylene components: **DO NOT USE IF THE STERILE PACKAGE APPEARS TO BE DAMAGED.**

For metal components only: If the sterile implant is determined to be aseptically compromised but still acceptable for intended use based on physician determination, the implant must be rinsed and sterilized prior to implantation according to the following instructions.

Rinsing/Cleaning

Use sterile room temperature water or physiological saline to soak the implant. Soak the implant for a minimum of 5 minutes. Immediately dry the product. Inspect the implant prior to sterilization.

Metal products that do not contain plastic or ceramic components and are not HA-coated or porous-coated may be resterilized, in a properly functioning, calibrated steam sterilizer. If sterilization of a **metal component** is necessary, the following parameters are recommended as they have been validated for a Sterility Assurance Level (SAL) of 10⁻⁶:

Method	Cycle	Temperature	Exposure Time
Steam	Pre-vacuum	270°F(132°C)	10 minutes

NOTE: The adequacy of any sterilization procedure must be suitably tested. It is critical that

appropriate process parameters be validated for each facility's sterilization equipment and product/load configuration by persons who have training and expertise in sterilization processes to substantiate the process and its reliability and reproducibility. Flash sterilization may be conducted, if applicable, according to the specific health care facility's policy.

The care and handling of porous-coated implants demands greater attention because of the increased potential for particulate and microbiological contamination. Body fluids, tissues, and particulate matter adhere to the beaded surfaces. Therefore, it is critical to minimize handling of the prosthesis.

The operating room staff should be conservative in component selection since opened packages will not be accepted for credit.

How Supplied

The NRH Knee and S-ROM Modular Total Knee Systems are supplied in a range of sizes. Please refer to the current price list for product codes and sizes available.

®S-ROM is a registered trademark of DePuy Orthopaedics, Inc.

™NOILES and NRH are trademarks of DePuy Orthopaedics, Inc.

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©S-ROM είναι σήμα κατατεθέν της DePuy Orthopaedics, Inc.
DePuy International Ltd.