**IMPORTANT INFORMATION – Please Read Before Use**  
**INFORMATIONS IMPORTANTES – À lire attentivement avant d'utiliser**  
**WICHTIGE INFORMATIONEN - Vor Gebrauch bitte lesen**  
**INFORMACIÓN IMPORTANTE - Sírvase leerla antes de emplear**  
**INFORMAZIONI IMPORTANTI - Leggere prima dell'uso**  
**INFORMAÇÃO IMPORTANTE - É favor ler antes de utilizar**  
**BELANGRIJKE INFORMATIE - Vóór gebruik lezen**  
**VIGTIG INFORMATION - Læses før brug**  
**VIKTIGT! - Läs detta innan produkten används**  
**TÄRKEÄÄ - Lue ennen käyttöä**  
**ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ - Παρακαλούμε διαβάστε πριν από τη χρήση**  
**DŮLEŽITÉ INFORMACE - Před použitím si prosím přečtete tento leták**  
**FONTOS INFORMÁCIÓ - Használat előtt elolvasandó**  
**WAŻNE INFORMACJE - Prosimy przeczytać przed użyciem**  
**DŮLEŽITÉ INFORMÁCIE - Prečítajte si láskavo ešte pred použitím**  
**TÄHTIS TEAVE – lugege enne kasutamist!**  
**SVARĪGA INFORMĀCIJA – Lūdzu, izlasiet pirms lietošanas**  
**SVARBI INFORMACIJA. Prašom perskaityti prieš naudojimą**  
**ÖNEMLİ BİLGİLER – Lütfen kullanmadan önce okuyun**  
**重要信息 – 请在使用前阅读**  
**重要資訊 – 使用前請詳細閱讀**  
**중요 정보 – 사용하기 전에 읽어 보시기 바랍니다!**  
**VIKTIG INFORMASJON – Må leses før bruk**

<b>Contents</b>	
English.....	5
<b>Contenus</b>	
Français.....	7
<b>Inhalt</b>	
Deutsch.....	9
<b>Contenido</b>	
Español.....	11
<b>Indice</b>	
Italiano.....	13
<b>Contém</b>	
Português.....	15
<b>Inhoud</b>	
Nederlands.....	17
<b>Indhold</b>	
Dansk.....	19
<b>Innehåll</b>	
Svenska.....	21
<b>Sisältö</b>	
Suomi.....	23
<b>Περιεχόμενα</b>	
Ελληνικά.....	25
<b>Obsah</b>	
Česky.....	27
<b>Tartalomjegyzék</b>	
Magyar.....	29
<b>Spis treści</b>	
Polski.....	31
<b>Obsah</b>	
Slovenčina.....	33
<b>Sisu</b>	
Eesti.....	35
<b>Saturs</b>	
Latviešu Valodā.....	37
<b>Turinys</b>	
Lietuvių Kalba.....	39
<b>İçindekiler</b>	
Türkçe.....	41
<b>目 录</b>	
.....	43
<b>目 録</b>	
.....	45
<b>목 차</b>	
.....	47
<b>Innhold</b>	
Norsk.....	49

INTERPRETATION OF SYMBOLS .....	51
SIGNIFICATION DES SYMBOLES.....	51
BEDEUTUNG DER SYMBOLE .....	51
INTERPRETACIÓN DE SÍMBOLOS.....	51
INTERPRETAZIONE DEI SIMBOLI.....	51
INTERPRETAÇÃO DE SÍMBOLOS.....	51
BETEKENIS VAN DE SYMBOLEN.....	51
BETYDNING AF SYMBOLER.....	51
FÖRKLARNING AV OLIKA SYMBOLERS BETYDELSE.....	51
SYMBOLIEN SELITYKSET.....	51
ΕΡΜΗΝΕΙΑ ΣΥΜΒΟΛΩΝ .....	51
VYSVĚTLIVKY VÝRAZŮ A SYMBOLŮ.....	51
MAGYARÁZAT A HASZNÁLT JELEKHEZ ÉS KIFEJEZÉSEKHEZ .....	51
ZNACZENIE SYMBOLI I TERMINÓW.....	51
VYSVETLENIE TERMÍNOV A SYMBOLOV .....	51
TINGMÄRKIDE TÄHENDUSED.....	51
SIMBOLU INTERPRETAČIJA.....	51
SIMBOLIŲ INTERPRETAVIMAS .....	51
SEMBOLLERİN ANLAMLARI.....	51
.....	51
.....	51
.....	51
SYMBOLENES BETYDNING.....	51

## ENGLISH

2

**STERILE**

See product label for method

### General Remarks

It is essential that the surgeon is fully conversant with the prescribed surgical technique and that due consideration is given to the information provided below. Improper patient selection, with particular regard to patient weight and functional demands, as well as placement and positioning of the implant may result in unusual stress conditions and subsequent reduction in their functional life.

Warning and instructions should be strictly observed.

Implants must not be modified or adapted unless they have been specifically designed to be so altered.

### Sterility

If DePuy implants are provided sterile, the integrity of the packaging should be checked to ensure that the sterility of the contents is not compromised.

If the packaging is damaged do not use or re-sterilise. Also unused components must not be re-sterilised.

### WARNINGS AND PRECAUTIONS

#### INFORMATION FOR USE

##### 1. Preoperative

The surgeon should discuss all physical and mental limitations particular to the patient and all aspects of the surgery and products with the patient before surgery. The discussion should include limitations of the procedure and the mechanical limitations of the material(s) selected for implantation. Factors that could impose limits on the performance and stability of the implant(s) e.g. activity level and patient weight, should be described to the patient so as to maximise the potential for long-term freedom from complications. The necessity of following the surgeon's instructions regarding these factors following the operation must be clearly understood by the patient.

An adequate inventory of sterile implant sizes should be on hand at the time of surgery.

##### 2. Handling

Extreme care must be taken while handling the implant so that it is not scratched or damaged in any way. Implants should be stored in their unopened original packaging. Where protective covers are provided they should be left in place until immediately prior to use. Implants are designed for single use and must never be re-implanted.

##### HA Coated Implants

HA coated implants must not be implanted with cement.

##### Modular Implants

DePuy modular implants featuring self-locking taper couplings are manufactured to extremely high tolerances. Under no circumstances should a DePuy modular implant component be combined with another component of a different manufacturer. Care must be taken to ensure that taper fit components are mutually compatible and that they are impacted together with the appropriate instruments. It is essential that compatible male and female tapers are examined to ensure they are completely clean prior to assembly.

##### Indications and Contraindications

Generally accepted indications for joint replacement include: severe arthropathy due to advanced rheumatoid or osteo-arthritis where conservative therapy or alternative treatments have failed or are considered unsuitable, arthropathy due to degenerative disease, acute trauma, and a previous failed joint replacement. The patient's joint must be anatomically and structurally suited to receive the selected implant(s).

Joint replacements may be contraindicated where the patient is overweight, where there is infection, poor bone stock, severe deformity, drug abuse, overactivity, tumour, mental incapacity, muscle, nerve or vascular disease.

##### 3. Surgical Technique

DePuy provides detailed written surgical instructions, videos and has access to other educational resources to ensure that the surgeon is conversant with the operative procedure. Further information on surgical techniques and DePuy products is available on request. Careful preoperative planning on the basis of radiographic findings should be carried out routinely. Radiographic templates are available for the majority of DePuy implants.

##### 4. Implant Fixation

When bone cement is used it is essential that the cement manufacturers' instructions for use are followed carefully. Handling of the bone cement can influence the efficacy of the implant fixation. For example, inclusion of blood and air and poor mixing of polymer and monomer may lead to early loosening or complications.

Care must be taken that the correct and appropriate size implant is used in conjunction with the correct instrumentation and trial components where these are provided. Definitive implants and trial components

manufactured by DePuy must not be used in conjunction with those of another manufacturer as component parts may not be compatible.

#### **5. Post-operative Care and Follow-up**

The patient should be released from the hospital with complete instructions and warnings, preferably in writing, regarding further exercises and therapies, and any limitations on their activities.

A continuing program of periodic follow-up controls is essential and should be strictly enforced. In every case accepted practices should be followed in post operative care.

The patient should be encouraged to report to his surgeon any unusual changes of the operated extremity. If evidence suggests loosening of the implant (particular pain and progressive changes in the radiographs) an intensified schedule of check-ups is advised and new warning and instructions to the patient may be appropriate regarding further activity restrictions.

In every case, accepted practices should be followed in post operative care. Excessive physical activity and trauma affecting the operated extremity have been implicated in the premature failure of joint arthroplasty as a result of position change, fracture and/or wear and tear of the implant. The functional life expectancy of prosthetic implants is at present not clearly established.

#### **Further information**

Further information is available from your DePuy representative on request.

#### **CAUTION**

Federal Law (USA) restricts this device to sale by or on the order of a physician.

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**D / DP** DePuy, Inc
**DePuy Orthopaedics, Inc.**
**DePuy International Ltd**