IMPORTANT INFORMATION – Please Read Before Use
INFORMATIONS IMPORTANTES - À lire attentivement avant d'utiliser
WICHTIGE INFORMATIONEN - Vor Gebrauch bitte lesen
INFORMACIÓN IMPORTANTE - Sírvase leerla antes de emplear
INFORMAZIONI IMPORTANTI - Leggere prima dell'uso
INFORMAÇÃO IMPORTANTE - É favor ler antes de utilizar
BELANGRIJKE INFORMATIE - Vóór gebruik lezen
VIGTIG INFORMATION - Læses før brug
VIKTIGT! - Läs detta innan produkten används
TÄRKEÄÄ - Lue ennen käyttöä
ΣΗΜΑΝΤΙΚΕΣ ΠΛΗΡΟΦΟΡΙΕΣ - Παρακαλούμε διαβάστε πριν από τη χρήση
DŮLEŽITÉ INFORMACE - Před použitím si prosím přečtete tento leták
FONTOS INFORMÁCIÓ - Használat előtt elolvasandó
WAŻNE INFORMACJE - Prosimy przeczytać przed użyciem
DÔLEŽITÉ INFORMÁCIE - Prečítajte si láskavo ešte pred použitím

Sterile	GP
----------------	-----------

 **DePuy**
a *Johnson & Johnson* company

0086

IFU-0902-00-710
Rev. B 12/04

**NEUFLEX™ MCP AND PIP FINGER JOINT IMPLANT
IMPLANTS NEUFLEX™ POUR ARTICULATIONS
METACARPO-PHALANGIENNES (MCP) ET
INTERPHALANGIENNES PROXIMALES (PIP)
NEUFLEX™ MCP- UND PIP-FINGERGELENK-IMPLANTATE
IMPLANTES DIGITALES NEUFLEX™ DE NUDILLOS MCP Y PIP
PROTESI PER L'ARTICOLAZIONE METACARPOFALANGEA (MCF) E
INTERFALANGEA PROSSIMALE (IFP) NEUFLEX™
IMPLANTE PARA AS ARTICULAÇÕES METACARPOFALÂNGICAS E
INTERFALÂNGICAS PROXIMAIS NEUFLEX™
NEUFLEX™ MCP- EN PIP-VINGERGEWRIGHT IMPLANTAATEN
NEUFLEX™ MCP OG PIP FINGERLEDSIMPLANTATER
NEUFLEX™ MCP- OCH PIP-FINGERLEDSIMPLANTAT
NEUFLEX™ MCP- JA PROKSIMAALINEN
INTERFALAGEAALINEN SORMINIVELIMPLANTTI
ΜΟΣΧΕΥΜΑΤΑ ΜΕΤΑΚΑΡΠΟΦΑΛΛΑΓΙΚΩΝ ΚΑΙ ΕΓΓΥΣ
ΦΑΛΑΓΓΟΦΑΛΑΓΓΙΚΩΝ ΔΑΚΤΥΛΙΚΩΝ ΑΡΘΡΩΣΕΩΝ NEUFLEX™
KLOUBNÍ IMPLANTÁT NEUFLEX™ PRO MCP A PIP SKLOUBENÍ PRSTŮ
NEUFLEX™ MCP ÉS PIP UJJÍZÜLETI IMPLANTÁTUM
PROTEZY STAWU PALCOWEGO NEUFLEX™ MCP ORAZ PIP FINGER
IMPLANTÁT METAKARPOFALANGÁLNEHO A PROXIMÁLNEHO
INTERFALANGÁLNEHO PRSTOVÉHO KLÍBA NEUFLEX™**

Contents		
English	Page	4
Contenus		
Français	Page	6
Inhalt		
Deutsch	Seite	8
Contenido		
Español	Página	11
Indice		
Italiano	Pagina	13
Índice		
Português	Página	16
Inhoud		
Nederlands	Pagina	19
Indhold		
Dansk	Side	21
Innehåll		
Svenska	Sida	23
Sisältö		
Suomi	Sivu	26
Περιεχόμενα		
Ελληνικά	Σελίδα	28
Obsah		
Česky	Str.	31
Tartalomjegyzék		
Magyar	Oldal	33
Spis treści		
Polski	Str.	36
Obsah		
Slovenčina	Str.	38
Interpretation of symbols	Page	41
Signification des symboles	Page	41
Bedeutung der Symbole	Seite	41
Interpretación de símbolos	Página	41
Interpretazione dei simboli	Pagina	41
Interpretação de símbolos	Página	41
Betekenis van de symbolen	Pagina	41
Betydning af symboler	Side	41
Förklaring av olika symbolers betydelse	Sida	41
Symbolien selitykset	Sivu	41
Ερμηνεία συμβόλων	Σελίδα	41
Vysvětlivky výrazů a symbolů	Str.	41
Magyarázat a használt jelekhez és kifejezésekhez	Oldal	41
Znaczenie symboli i terminów	Str.	41
Vysvetlenie termínov a symbolov	Str.	41

**DESCRIPTION**

The NeuFlex™ MCP and PIP Finger Joint Prostheses are flexible, one-piece, hinged implants manufactured from silicone elastomer. The NeuFlex MCP is designed to replace a damaged metacarpophalangeal (MCP) joint and the NeuFlex PIP is designed to replace a damaged proximal interphalangeal (PIP) joint. The proximal and distal stems of the prosthesis form an angle which mimics the approximate position of the joint when the hand is relaxed.

INDICATIONS

The NeuFlex MCP and NeuFlex PIP Finger Joint Prostheses are indicated for cementless replacement of the metacarpophalangeal (MCP) and interphalangeal (PIP) joints, respectively, where disabled by rheumatoid, degenerative or traumatic arthritis.

CONTRAINDICATIONS

The following conditions are relative contraindications for use of this prosthesis:

1. Active local or systemic infection;
2. Destruction of the metacarpal or phalanx or poor bone quality which prevents adequate fixation of the implant;
3. Loss of musculature, neuromuscular compromise or vascular deficiency in the affected finger;
4. Growing patients with open epiphyses;
5. Patients with high activity levels;
6. Patients unwilling or unable to comply with the physician's instructions.

STERILITY

The NeuFlex Finger Joint Prosthesis is supplied **sterile**. Remove from the package using accepted aseptic technique, only after the correct size has been determined.

This product should not be used if the sterile packaging is damaged such that sterility may be compromised.

Resterilization of the prosthesis is not recommended.

HANDLING

Proper handling of any implant is mandatory. To avoid nicking, tearing, cutting or scratching this silicone rubber, it should be handled only with blunt instruments and should never be handled with instruments with sharp edges, teeth or serrations. Prior to surgical use, each implant should be inspected visually for possible imperfections. Once removed, an implant should never be re-implanted. Even if the implant appears undamaged it may be fatigued from previous use and may have developed microscopic imperfections which could lead to implant failure.

PREOPERATIVE

THE SURGEON SHOULD DISCUSS ALL PHYSICAL AND MENTAL LIMITATIONS PARTICULAR TO THE PATIENT AND ALL ASPECTS OF THE SURGERY AND PRODUCT WITH THE PATIENT BEFORE SURGERY. The discussion should include the limitations of joint reconstruction, the possible consequences resulting from these limitations and the necessity of following the surgeon's instructions postoperatively, particularly with regard to rehabilitation and patient activity levels. Multiple sizes of sterilized implants should be on hand at the time of surgery. It is recommended that this include at least one size larger and

one size smaller than the preoperatively determined size to accommodate intraoperative selection of an appropriate size implant.

INTRAOPERATIVE

Trial prostheses should be used to determine the correct size of implant. The implant itself should not be used as a trial implant.

There is evidence that the potential for sepsis following joint arthroplasty may be reduced by:

1. Consistent use of prophylactic antibiotics;
2. Having all operating room personnel (anesthetists, circulators, observers, etc.), as well as scrubbed technicians, properly attired; and
3. Utilization of a laminar flow clean air system within the operating room.

Familiarity with and careful attention to the surgical technique are necessary for good results. Improper selection, placement, or positioning of the implant may result in unusual stress conditions and subsequent reduction in the functional life. **IMPLANTS AND TRIAL IMPLANTS FROM ONE MANUFACTURER SHOULD NOT BE USED TOGETHER WITH THOSE OF ANOTHER MANUFACTURER SINCE MATCHING OF DIFFERENT COMPONENTS CANNOT BE ASSURED.** Prior to closure, the surgical site should be thoroughly cleaned and range of motion should be thoroughly checked for implant impingement or instability.

POSTOPERATIVE

For the postoperative management and gradual functional recovery of the operated extremity, strict adherence by the patient to the surgeon's instructions and warnings is extremely important. The operated extremity should be stabilized appropriately. The patient should be released from the hospital with complete instructions and warnings, preferably in writing, regarding rehabilitation exercises and therapies and any limitations on their activities. A continuing program of periodic patient follow-up is recommended. In every case, accepted practices should be followed in postoperative care. Excessive physical activity and trauma affecting the replaced joint have been implicated in premature failure of the reconstruction by position change, fatigue failure, and/or wear of the implant. The functional life expectancy of silicone elastomer implants is at present not clearly established.

ADVERSE EFFECTS AND COMPLICATIONS

The following are generally the most frequent adverse effects or complications encountered in finger joint replacement. These adverse effects can often be related to factors listed under **"WARNINGS AND PRECAUTIONS"**.

1. Failure of the implant due to fatigue, wear or over-loading;
2. Early or late infection;
3. Wear particles caused by the movement and wear of a silicone rubber implant may cause or exacerbate synovitis or bone cyst formation.
4. There have been reports in the literature which suggest that some individuals may have an immunological reaction to silicone implants, resulting in connective tissue and/or autoimmune disease. If these conditions are suspected, removal of this silicone device should be considered.

WARNINGS AND PRECAUTIONS

The following conditions, singularly or concurrently, tend to place excessive loads on a finger joint implant and, thereby, place the patient at higher risk for failure of finger joint replacement. **If excessive loading of the affected finger joint cannot be prevented, a finger joint implant should not be used.**

1. excessive activity of the affected joint;
2. uncorrected or recurrent deformity;

3. incorrect sizing of the implant;
4. inadequate soft tissue or bony support;
5. implant malposition.

The benefits of finger joint replacement may not meet the patient's expectations or may deteriorate over time. Pain, swelling, instability and/or deformity may persist or return after finger joint replacement.

INFORMATION

Contact your DePuy representative regarding surgical technique brochures or for further information about this implant.

CAUTION

Federal Law (USA) restricts this device to sale by or on the order of a physician.

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DePuy Orthopaedics, Inc.
DePuy International Ltd.