

(GB) Avicenum PHLEBO 360 medical compression products (CCL 2)

(DE) Avicenum PHLEBO 360 Medizinische Kompressionsprodukte (CCL 2)

(RU) Изделия медицинские компрессионные антиварикозные Avicenum 360

(PL) Medyczne kompresyjne produkty Avicenum PHLEBO 360 (CCL 2)

(CZ) Zdravotní kompresivní výrobky Avicenum PHLEBO 360 (CCL 2)

(SK) Zdravotné kompresívne výrobky Avicenum PHLEBO 360 (CCL 2)



LYCRA®Soft Comfort




Covered LYCRA®



Microfibre



EDEMA ELIMINATION  
ZUR VERRINGERUNG VON ÖDEMEN  
УСТРАНЕНИЕ ОТЕКОВ  
ZŁAGODZENIE OBRZĘKÓW  
ELIMINÁCE OTOKU  
ELIMINÁCIA OPŮCHOV



GRADUATED COMPRESSION  
MIT KOMPRESSIONSGRADIENT  
ГРАДИЕНТ СЖАТИЯ  
STOPNIOWA KOMPRESJA  
STUPŇOVANÁ KOMPRESIE  
STUPŇOVANÁ KOMPRESIA



IMPROVEMENT OF VENOUS RETURN  
VERBESSERTER VENÖSER BLUTFLUSS  
УЛУЧШЕНИЕ ВЕНОЗНОГО ВОЗВРАТА  
POPRAWA POWROTU ŻYLNĘGO  
ZLEPŠENÍ ŽILNÍ NÁVRATNOSTI  
ZLEPŠENIE ŽILOVEJ NÁVRATNOSTI



SUPPORT OF BLOOD CIRCULATION  
BLUTKREISLAUF UNTERSTÜTZUNG  
ПОДДЕРЖКА КРОВООБРАЩЕНИЯ  
WSPARCIE KRWIOBIEGU  
PODPORA KREVŇNÍHO OBĚHU  
PODPORA KRVNÉHO OBĚHU



**Aries**  
MEDICAL PRODUCTS

www.aries.eu



ARIES, a.s.  
Studenec 309, 512 33  
MADE IN CZECH REPUBLIC  
ISO 9001:2015



*Avicenum*  
360



... for easier life!



Medical  
compression products



Medizinische  
Kompressionsprodukte



Изделия  
медицинские  
компрессионные  
антиварикозные



Medyczne  
kompresyjne produkty



Zdravotní  
kompresivní výrobky



Zdravotné  
kompresívne výrobky

**Aries**  
MEDICAL PRODUCTS

*Avicenum*  
360



**PHLEBO**

Medical  
compression  
products

**FINE**

CCL  
2

23-32 mmHg



Sanitized®

Deo and bacteriostatic effect.

*Avicenum*  
360



... for easier life!

**Material:**

60 % PAD Nylon/PAD  
Nylon/PAD Нейлон/PAD  
Nylon/PAD Nylon/PAD Nylon  
40 % elastane LYCRA®/  
Elasthan LYCRA®/эластан  
LYCRA®/LYCRA® dzianina  
elastyczna/elastan  
LYCRA®/elastan LYCRA®

**Sanitized® composition:**

Dimethyltetradecyl  
(3-(trimethoxysilyl)propyl)  
ammonium chloride

Do not use fabric  
softeners/Keinen Weich-  
püler verwenden/Не исполь-  
зовать кондиционер для белья/  
Nie używać płynów zmiękc-  
zających/Nepoužívať aviváž/  
Nepoužívať aviváz



13916\_05\_2021