

Vyobrazenie zdravotníckej pomôcky

Systém koronárneho bifurkačného stentu Biomime Branch uvoľňujúci sirolimus - rôzne veľkosti

P6850A

biomime™ branch
Sirolimus Eluting Coronary Side Branch Stent System
Intuitive Solution for Bifurcations

BioMime Branch Stent Specifications

| Stent Length (LL) | Length of L1 | Length of L2 | Length of L3 |
|-------------------|--------------|--------------|--------------|
| 16 mm | 5.0 mm | 4.5 mm | 6.5 mm |
| 19 mm | 5.5 mm | 4.5 mm | 9.0 mm |
| 24 mm | 6.5 mm | 4.5 mm | 13.0 mm |
| 29 mm | 6.5 mm | 4.5 mm | 18.0 mm |

BioMime Branch Deployment Technique

STEP 1
Prepare and Predilate the side branch lesion

STEP 2
Deployment of BioMime Branch

STEP 3
Post dilation

STEP 4
Prepare the main branch

STEP 5
Deploy stent in main branch and follow up with Kissing Balloon Technique with both main and side branch