

May 2023



Important Update on ZEISS toric IOLs – Conversion to Spherical Equivalent Notation

Dear customer,







We are pleased to inform you about the upcoming switch of the notation from *Sphere (SPH)* to *Spherical Equivalent (SE)* in the ZEISS toric IOL portfolio.

The **AT TORBI 709M/MP** and **AT LISA tri toric 939M/MP** are labelled in SPH/Cyl and will be switched to the SE/Cyl models **AT TORBI 719M/MP** and **AT LISA tri toric 949M/MP**.

Why are we doing this?

After switching the notation from Sphere (SPH) to Spherical Equivalent (SE) the ZEISS toric IOL portfolio will be fully streamlined with the ZEISS CATARACT Workflow, thereby **simplifying ZEISS IOL Calculation and ZEISS IOL Ordering**.

What's new?

	Before Switch	After Switch
Product Name	AT TORBI® 709M/MP and AT LISA® tri toric 939M/MP	AT TORBI 719M/MP and AT LISA tri toric 949M/MP
Product Label	 SPH and CYL notation	 SE / CYL notation only Simplified design of labels with data matrix code
Package Design		 New look and easier lens type identification
Product Feature	 <ul style="list-style-type: none">AT TORBI 709M/MP■ No anti-PCO ring■ No orientation holes■ Company labeledA-Constant: 118.3	 <ul style="list-style-type: none">AT TORBI 719M/MP■ 360 degree anti-PCO ring around the posterior optic■ Orientation holes■ Company labeledA-Constant: 118.4

Sitz der Gesellschaft:
Goeschwitzer Strasse 51 - 52
07745 Jena

Lieferanschrift:
Carl Zeiss Meditec AG
Max-Dohrn-Strasse 8 - 10
10589 Berlin

Banken:
Deutsche Bank Jena
Konto: 624536900 (BLZ 820 700 00)
IBAN: DE90 8207 0000 0624 5369 00
BIC/ SWIFT: DEUT DE 8EXXX

Commerzbank Jena
Konto: 258072800 (BLZ 820 400 00)
IBAN: DE31 8204 0000 0258 0728 00
BIC/ SWIFT: COBADEFFXXX

Amtsgericht Jena HRB 205623
UST-ID. Nr.: DE 811 922 737
WEEE-Reg.-Nr.: DE55298748

Vorsitzender des Aufsichtsrats:
Dr. Karl Lamprecht

Der Vorstand:
Dr. Markus Weber (Vorsitzender)
Justus Felix Wehmer

To ensure a smooth transition, ZEISS is supporting the switch as following:

- **Automated update:** IOL Calculation Service and ZCALC will calculate only the SE models on the switch date.
- **Assisted update:** IOLMaster 500 / 700 and EQ Workplace - The ZEISS IOL sales team will support you directly during the required updates.


Please note that the SPH models cannot be ordered anymore after the switch date.

If you have questions or need further assistance, please do not hesitate to reach out to your ZEISS IOL sales representative.


Best regards,

Your ZEISS Toric IOL Team

Carl Zeiss Meditec AG



i.A. Bernd Kimmerle
Head of Product Management Refractive IOLs



i.A. Kim, En Gyung
Product Manager Refractive IOLs