

Angiographic Catheter

- Cathéter pour angiographie ■ Catéter para angiografia ■
- Angiographiekatheter ■ Catetere per angiografia ■
- Cateter angiográfico ■ Angiografikater ■
- Angiografikateter ■ Angiografia Catheter ■
- Αγγιογραφικός Καθετήρας ■ 血管造影用カテーテル ■


Accu-VU™

NON-BRAIDED

Sans treillis / Sin malla interna / Ohne Gewebeverstärkung /
Non intrecciato / Sem tranção / Non-Braided /
Ikke-Flettet / Χρόμη Πλέγμα / 非編組み

Sizing Pigtail

with (2) R/O Markers
(2cm)



SIZE / Taille / Variante / Größe / Misure /
Tamanho / Storlek / Størrelse / Maat / Μέγεθος / サイズ

4 F

LENGTH / Longueur / Longitud / Länge /
Lunghezza / Comprimento / Längd / Længde /
Lengto / Μήκος / 長さ

70 cm

DESIGNED TO ACCEPT
GUIDEWIRE / Conçu pour recevoir un guide /
Guia recomendada / Akzeptiert einen Führungsdraht /
Da usare con guida / Projelado para receber fio-guida /
Anpassad för ledare / Tilpasset til leder / Passende Voerdraad /
Συμβατικό με οδηγό / Οδηγός / 適合ガイドワイヤーに適合

.035 in.

SIDEHOLES / Ouvertures latérales /
Orificios laterales / Seitenlöcher / Fori laterali /
Orificios laterais / Sidohal / Sidehuller / Zijgaten /
Παραπλευρικές οπές / サイドホール

10

MAXIMUM FLOW RATE

Débit maximal / Flujo máximo / Maximal
Flußgeschwindigkeit / Anjenida de máxima
de fluxo / Ritmo de fluxo máximo
Maximal flödeshastighet / Maksimal
flöshastighed / Μέγιστος ρυθμός ροής
最大流量

20cc/sec

PRESSURE LIMIT

Pression maximale / Limite de pression
Maximaler Druck / Pressione massima
Limite De Pressão / Trykgrænse
/ Trykgrænse / Μέγιστο ΠΕΑ / 耐圧限度

1050psi

Catalog No. REF

13709105

LOT

521057

Use By


2013-08

QTY

1

Quantité / Cantidad / Menge /
Quantia / Quantidade / Antal / Antal /
Hoeveelheid / Ποσότητα / 数量

- DO NOT REUSE • STERILIZED USING ETHYLENE OXIDE • DO NOT RESTERILIZE
- CAUTION • KEEP DRY • UPPER TEMPERATURE LIMIT • KEEP AWAY FROM SUNLIGHT
- DO NOT USE IF PACKAGE IS DAMAGED • STERILE UNLESS PACKAGE IS OPENED OR DAMAGED



CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

MAXIMUM FLOW RATE

Débit maximal / Flujo máximo / Maximal
Flußgeschwindigkeit / Anjenida de máxima
de fluxo / Ritmo de fluxo máximo
Maximal flödeshastighet / Maksimal
flöshastighed / Μέγιστος ρυθμός ροής
最大流量

20cc/sec

PRESSURE LIMIT

Pression maximale / Limite de pression
Maximaler Druck / Pressione massima
Limite De Pressão / Trykgrænse
/ Trykgrænse / Μέγιστο ΠΕΑ / 耐圧限度

1050psi

Catalog No. REF

13709105

LOT

521057

Use By


2013-08

QTY

1

Quantité / Cantidad / Menge /
Quantia / Quantidade / Antal / Antal /
Hoeveelheid / Ποσότητα / 数量

- DO NOT REUSE • STERILIZED USING ETHYLENE OXIDE • DO NOT RESTERILIZE
- CAUTION • KEEP DRY • UPPER TEMPERATURE LIMIT • KEEP AWAY FROM SUNLIGHT
- DO NOT USE IF PACKAGE IS DAMAGED • STERILE UNLESS PACKAGE IS OPENED OR DAMAGED



CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician.

CE
0086

EC REP

AngioDynamics UK Ltd Building 2000,
Beach Drive IQ Cambridge, Waterbeach,
Cambridge, CB25 9TE United Kingdom
Telephone: +44 (0) 1223 729300
Fax: +44 (0) 1223 729329

MANUFACTURED IN USA ANGIODYNAMICS, INC. LC335 REV.03

Rev. A
MADE IN THE USA
U.S. PATENT: 5,584,821

AngioDynamics, Inc.
603 Queensbury Avenue
Queensbury, NY USA 12804

For orders and customer service in the U.S.A. call
toll free: 1-800-772-6446 or fax: 1-518-798-1360



ANGIODYNAMICS INC CATALOG No 13709105

007