


The following is a representative sample of the final printed label. Verify placement of information according to this sample. DO NOT refer to this sample for LOT or DATE information.

Armada 14 Label



# Armada 14

## PTA Catheter

**OTW**

**.014**

**Armada 14**


**1.5 20 90**  
mm mm cm

**Armada 14**  
PTA Catheter

REF A1015-020

1.5 mm x 20 mm x 90 cm

**BALLOON** **BALLOON** **USAB CAT**



REC MIN INTRO 4 F ID 0.053" 1.35 mm

REC MIN INTRO 4 F ID 0.053" 1.35 mm

85°F 30°C

15°C

WPL2093904-01 (4/2/13)

REF A1015-020 2016-09 LOT 30411T1

(01)09717648155628(17)160930(10)30411T1(91)3000

**REC MIN INTRO**

**REC GW**

**NPR**

**RBP**

**4 F**  
ID 0.053"  
1.35 mm

**0.014"**  
0.35 mm

**8 atm**  
800 kPa

**14 atm**  
1400 kPa

| Balloon Compliance |         |      |
|--------------------|---------|------|
| Pressure           | BALLOON |      |
| [atm]              | [kPa]   | [mm] |
| 4                  | 400     | 1.48 |
| 5                  | 500     | 1.47 |
| 6                  | 600     | 1.48 |
| 7                  | 700     | 1.49 |
| *8                 | 800     | 1.50 |
| 9                  | 900     | 1.51 |
| 10                 | 1000    | 1.52 |
| 11                 | 1100    | 1.53 |
| 12                 | 1200    | 1.54 |
| 13                 | 1300    | 1.55 |
| **14               | 1400    | 1.56 |

\* NPR

\*\* RBP DO NOT EXCEED!

**Rx ONLY** **CE** 0123

**STERILE EO**

**NON PYROGENIC**

2013-04

Legal Manufacturer  
Abbott Laboratories Vascular Enterprises  
Anthoptstrasse  
8222 Beringen - Switzerland  
Phone: +1 800 227-9902  
FAX: +1 800 601-8874  
Outside USA Phone: +1 951 914-4888  
Outside USA FAX: +1 951 914-2531

**MANUFACTURED BY**  
Abbott Vascular  
Cashel Road  
Clonmel  
County Tipperary, Ireland

Abbott Laboratories Vascular Enterprises  
Armada 14  
1.5 mm x 20 mm x 90 cm  
REF A1015-020  
LOT 30411T1 OTW

Abbott Laboratories Vascular Enterprises  
Armada 14  
1.5 mm x 20 mm x 90 cm  
REF A1015-020  
LOT 30411T1 OTW

**Armada 14**  
**BALLOON** 1.5 mm **BALLOON** 20 mm **USAB CAT** 90 cm

(01) 0 8717648 15562 8

(17) 160930 (10) 30411T1 (91) 3000

LOT 30411T1 2013-04 2016-09 OTW

RELEASED

CONFIDENTIAL


May not be reproduced outside Abbott Vascular without written permission from Document Control.

The following is a representative sample of the final printed label. Verify placement of information according to this sample. DO NOT refer to this sample for LOT or DATE information.

**Armada 35 Label**

# Armada 35

## PTA Catheter



OTW

OTW

OTW

OTW

**.035**

Armada 35

3.0

20

80

BALLOON

BALLOON

USAB CAT

REF B1030-020

2016-07

LOT 30213T1

### Armada 35

#### PTA Catheter


REF B1030-020

3.0 mm x 20 mm x 80 cm

BALLOON

BALLOON

USAB CAT





2013-02

REC MIN INTRO

REC GW

NPR

RBP

6 Fr  
ID 1.78 mm

0.036"   
0.89 mm

8 atm  
800 kPa

28 atm  
2800 kPa

R&CE ONLY

0123

STERILE R

NON PYROGENIC

Abbott Laboratories  
Vascular Enterprises  
Armada 35  
3.0 mm x 20 mm  
x 80 cm  
REF B1030-020  
LOT 30213T1 OTW

Abbott Laboratories  
Vascular Enterprises  
Armada 35  
3.0 mm x 20 mm  
x 80 cm  
REF B1030-020  
LOT 30213T1 OTW

**Armada 35**  
BALLOON 3.0 mm    BALLOON 20 mm    USAB CAT 80 cm  
(01) 0 8717648 15448 5  
(17) 160731 (10) 30213T1 (91) 3000  
LOT 30213T1    2013-02    2016-07    OTW

Legal Manufacturer  
Abbott Laboratories Vascular Enterprises  
Anthoptstrasse  
8222 Beringen, Switzerland  
Phone: +1 800 227-9802  
FAX: +1 800 801-8874  
Outside USA Phone: +1 951 914-4669  
Outside USA FAX: +1 951 914-2531  
  
**MANUFACTURED BY**  
Abbott Vascular  
Coyol Free Zone  
El Coyol Alajuela,  
Costa Rica

011671-001/004

**CONFIDENTIAL**

May not be reproduced outside Abbott Vascular without written permission from Document Control.



Armada 35 LL Label

OTW

OTW

OTW

OTW

# Armada 35 LL

## PTA Catheter

.035

Armada 35 LL

Armada 35 LL

PTA Catheter

REF B1040-150

4.0 mm x 150 mm x 80 cm

BALLOON
BALLOON
USAB CAT

4.0
150
80

BALLOON mm
BALLOON mm
USAB CAT cm

REC MIN INTRO

REC GW

NPR

RBP

5 Fr  
ID 1.78 mm

0.035"  
0.89 mm

8 atm  
800 kPa

20 atm  
2000 kPa

REF B1040-150

2016-07

LOT 30213T1

(01)08717648154560(17)160731(10)30213T1(91)3000

REC MIN INTRO 16 Fr ID 1.78 mm

2013-02

Abbott Laboratories Vascular Enterprises

Armada 35 LL

4.0 mm x 150 mm x 80 cm

REF B1040-150

LOT 30213T1

Abbott Laboratories Vascular Enterprises

Armada 35 LL

4.0 mm x 150 mm x 80 cm

REF B1040-150

LOT 30213T1

WPL2088825-09 (2/11/13)

Legal Manufacturer: Abbott Laboratories Vascular Enterprises

Anthoptstrasse 8222 Beringen, Switzerland

Phone: +1 800 227-9902

FAX: +1 800 601-9874

Outside USA Phone: +1 951 914-4669

Outside USA FAX: +1 951 914-2531

CONFIDENTIAL

May not be reproduced outside Abbott Vascular without written permission from Document Control.

RELEASED - EFFECTIVE: 24 Jan 2014



# Armada 14 XT

## PTA Catheter

OTW

.014

**Armada 14 XT**

**1.20 8 145**  
mm mm cm

REF F2012-008 LOT 31202T1 2016-11




(01)08717648182808(17)161130(10)31202T1(91)3000

### Armada 14 XT

#### REF F2012-008

1.20 mm x 8 mm x 145 cm

PTA Catheter; PTA-Katheter; Cathéter d'ATP; Catéter de ATP; Catetere per PTA; Cateter para ATP; PTA-kateter; PTA-katheter; PTA-kateter; Καθετήρας PTA; PTAカテーテル; Cewnik do PTA; PTA katéter; Katetr pro PTA; PTA Kateteri; PTA-kateter; Катетер за PTA; Cateter pentru ATP; Катетер для ЧТА; Pallolaajennuskatetri; Katéter PTA; PTA kateeter; PTA katetrs; PTA kateteris; kateter za PTA



Consult instructions for use; Gebrauchsanweisung lesen; Consulter le mode d'emploi; Consultar las instrucciones de uso; Consultare le istruzioni per l'uso; Consultar as instruções de utilização; Se bruksanvisningen; Raadpleeg gebruiksaanwijzing; Lees bruksanvisningen; Συμβουλευτείτε τις οδηγίες χρήσης; 添付文書をご参照下さい; Zapoznać się z instrukcją obsługi; Lásd a használati utasítást; Viz návod k použití; Kullanım talimatlarına bakın; Les bruksanvisningen; Вижте инструкциите за употреба; Consultați instrucțiunile de utilizare; Ознакомьтесь с инструкциями по применению; Katso käyttöohjeita; Pozri návod na použitie; Vaadake teavet kasutusjuhistest; Skatiet lietošanas instrukciju; Žr. naudojimo instrukcijas; Glejte navodila za uporabo



**STERILE EO**



WPL2091586-02 (11/18/13)

Armada is a trademark of the Abbott Group of Companies.

| ATM      | kPa  | ↔       |
|----------|------|---------|
| 2        | 200  | 1.10 mm |
| 3        | 300  | 1.12 mm |
| 4        | 400  | 1.13 mm |
| 5        | 500  | 1.14 mm |
| 6        | 600  | 1.17 mm |
| 7        | 700  | 1.18 mm |
| 8 (NOM)  | 800  | 1.20 mm |
| 9        | 900  | 1.21 mm |
| 10       | 1000 | 1.23 mm |
| 11       | 1100 | 1.24 mm |
| 12       | 1200 | 1.25 mm |
| 13       | 1300 | 1.26 mm |
| 14 (RBP) | 1400 | 1.27 mm |
| 15       | 1500 | 1.28 mm |
| 16       | 1600 | 1.29 mm |
| 17       | 1700 | 1.30 mm |

**Legal Manufacturer**  
Abbott Vascular  
3200 Lakeside Drive  
Santa Clara, CA 95054 USA  
TEL: (800) 227-9902  
FAX: (800) 601-8874  
Outside USA TEL: (951) 914-4669  
Outside USA FAX: (951) 914-2531

#### MANUFACTURED BY

Abbott Vascular  
Coyol Free Zone  
El Coyol Alajuela, Costa Rica

#### EC REP

Abbott Vascular International BVBA  
Park Lane, Culliganlaan 2B  
1831 Diegem, BELGIUM  
TEL: + 32 2 714 14 11  
FAX: + 32 2 714 14 12



PPL2083844 (12/12/11)

Abbott Vascular  
Armada 14 XT  
1.20 mm x 8 mm  
x 145 cm  
REF F2012-008  
LOT 31202T1 OTW

Abbott Vascular  
Armada 14 XT  
1.20 mm x 8 mm  
x 145 cm  
REF F2012-008  
LOT 31202T1 OTW

Abbott Vascular  
Armada 14 XT



(01) 0 8717648 18280 8

REF F2012-008



(17) 161130 (10) 31202T1 (91) 3000

LOT 31202T1



2013-12



2016-11

Date printed: Máj 06, 2014

VERIFY THAT THIS IS THE CORRECT REVISION BEFORE USE