

Cera™

1012

Patent Ductus Arteriosus Occluder

EN Patent Ductus Arteriosus Occluder	FR Fermeture Persistence du Canal Artériel (PCA)
DE Patentductus Arteriosus-Verschlussystem	
IT Pervietà Dotto Arterioso Occluder	ES Ocluser para comunicaci3n Arterial permanente
PT Dispositivo de encerramento de comunica33es arteriais patentes (PDA)	SV Patent Ductus Arteriosus Occluder
TR Patent Duktus Arteriyozus Kapatma Cihazı	NL Patent Ductus Arteriosus Sluiting Systeem


REF LT-PDA-1012 LOT 10121611

2010-12 2013-12 SN 001

Proximal Dia.(ØD1) 10 mm

Distal Dia.(ØD2) 12 mm

Waist Length(L) 7 mm



Patent Ductus Arteriosus Occluder

REF 2013-12

REF LT-PDA-112

Cera™



Cera™ PDA Occluder



REF LT-PDA-1012

2013-12



Cera™ PDA Occluder



REF LT-PDA-1012

2013-12



- Consult instructions for use
- Do not reuse
- Store at room temperature
- Sterilized using ethylene oxide
- Non-pyrogenic
- Keep dry
- Do not use if package is damaged
- Do not resterilize
- Keep away from sunlight

Lifetech Scientific (Shenzhen) Co., Ltd.
 Floor 1-3, Cybio Electronic Building, Langshan 2nd Street,
 North Area of High-tech Park, Nanshan District,
 Shenzhen 518057, P.R.China
 Tel:+86-755-86026250 Fax:+86-755-86026251

Wellkang Ltd (www.CE-marking.eu)
 29 Harley St., London W1G 9QR, UK
 Tel:+44(20)88168309 Fax:+44(20)76811874



0344

V1.3 04/2011

Cera™

1012

Patent Ductus Arteriosus Occluder

EN Patent Ductus Arteriosus Occluder	FR Fermeture Persistence du Canal Artériel (PCA)
DE Patentductus Arteriosus-Verschlussystem	
IT Pervietà Dotto Arterioso Occluder	ES Ocluser para comunicaci3n Arterial permanente
PT Dispositivo de encerramento de comunica33es arteriais patentes (PDA)	SV Patent Ductus Arteriosus Occluder
TR Patent Duktus Arteriyozus Kapatma Cihazı	NL Patent Ductus Arteriosus Sluiting Systeem

REF LT-PDA-1012 LOT 10121611

2010-12 2013-12 SN 001



Cera™ PDA Occluder



REF LT-PDA-1012

2013-12



- Consult instructions for use
- Do not reuse
- Store at room temperature
- Sterilized using ethylene oxide
- Non-pyrogenic
- Keep dry
- Do not use if package is damaged
- Do not resterilize
- Keep away from sunlight



0344

Lifetech Scientific (Shenzhen) Co., Ltd.
 Floor 1-3, Cybio Electronic Building, Langshan 2nd Street, North Area of
 High-tech Park, Nanshan District, Shenzhen 518057, P.R.China
 Tel:+86-755-86026250 Fax:+86-755-86026251

Wellkang Ltd (www.CE-marking.eu)
 29 Harley St., London W1G 9QR, UK
 Tel:+44(20)88168309 Fax:+44(20)76811874

V1.3 04/2011