

Vyobrazenie zdravotníckej pomôcky

Katéter RX PTCA balónikový dilatačný uvoľňujúci sirolimus Mozec SEB -

rôzne veľkosti

P6851A



MOZEC™ SEB
Sirolimus Eluting Rx PTCA Balloon Dilatation Catheter

Handicraft for *Angiocraft*.

Nanotech Novel Drug Delivery

- Unique formulation - Solid lipid nano spheres (SLN) consisting of Sirolimus + lipid
- Drug concentration Sirolimus 3.00 µg/mm²
- Stable drug formulation with controlled & targeted drug release
- Rapid release of drug into tissue with drug reservoir effect, different size nano spheres allow drug delivery to each layer

Deliverability enhanced with Mozec PTCA catheter

- Low tip entry profile and semi-compliant balloon material allowing flat compliance
- PTFE coated proximal shaft for enhanced navigability and to cross challenging lesions
- Transparent ergonomic hub displaying balloon dimensions with flexible distal segment and pushable proximal segment

Size Matrix

Ø - 2.00, 2.25, 2.50, 2.75, 3.00, 3.50, 4.00, 4.50 mm
/ - 9, 12*, 14, 15*, 17, 20, 25, 30 mm

* For 200 mm balloon diameter only