


Effecta D

Dual Chamber Pacemaker

- Atrial & Ventricular Capture Control
- IRS^{plus}
- AutoSensing®
- Auto-Initialization
- Follow-up Center



Housing

Dimensions/weight 53 x 43 x 6.5 mm/26 g

Volume 11 cm³

Electrically conductive housing surfaces

➤ Uncoated 33 cm²

➤ Coated 7 cm²

X-ray identification SF

Ordering information

➤ Uncoated 375 429

➤ Coated 375 428

- 1) Only available for mode switching
- 2) EN 50061 triangle pulse
- 3) If Capture Control is ON, the pulse amplitude is automatically selected

- 4) 300...[25]...**350**...[25]...775 ms for AAI, AAT, DDT modes
- 5) Post-ventricular atrial blanking
- 6) Storage of IEGMs by using intelligent memory management
- 7) See manual for other modes
- 8) Nominal data of the manufacturer

All data at 37°C, 500 Ω.
Default settings are printed in bold.

BIOTRONIK SE & Co. KG
 Woermannkehre 1
 12359 Berlin · Germany
 Tel +49 (0) 30 68905-0
 Fax +49 (0) 30 6852804
 sales@biotronik.com
 www.biotronik.com