



Cardiac Rhythm Management

Tachyarrhythmia Therapy

Single-chamber ICD

Ilesto 5 VR-T

MR Conditional single-chamber ICD

ProMRI®



Product Highlights

■ **ProMRI®**

Allows patients to undergo MR scanning under specific conditions.

■ **Small size**

Increases the patients' comfort through a reduced device thickness.

■ **SafeSync RF telemetry**

RF telemetry for wandless, time-saving and reliable data transmission at implantation and follow-up.

■ **Ventricular Capture Control**

Automatic adjustment of pacing amplitudes for improved patient safety and extended device longevity.

■ **SelectSense Advanced**

Adaptation of sensing parameters to patients' individual needs via a sophisticated automatic sensitivity control (ASC) algorithm and several preset options.

■ **DF4 connectors**

For simplified and shortened implantation procedures and less material in the device pocket.

■ **BIOTRONIK Home Monitoring®**

Enables wireless patient monitoring, including Home Monitoring-supported follow-ups with remote scheduling, for early detection of clinical and device-related events.

Ordering information

Model	Connectors	Volume/weight	Dimensions	Order number
Ilesto 5 VR-T	DF-1 (2x) IS-1 (1x)	33 cm ³ /80 g	65 mm x 55 mm x 11 mm	390118
Ilesto 5 VR-T	DF4 (1x)	31 cm ³ /80 g	65 mm x 52 mm x 11 mm	390120



Cardiac Rhythm Management

Tachyarrhythmia Therapy

Single-chamber ICD

Ilesto 7 VR-T

MR Conditional single-chamber ICD

ProMRI®



Product Highlights

■ **ProMRI®**

Allows patients to undergo MR scanning under specific conditions.

■ **Small size**

Increases the patients' comfort through a reduced device thickness.

■ **SafeSync RF telemetry**

RF telemetry for wandless, time-saving and reliable data transmission at implantation and follow-up.

■ **Ventricular Capture Control**

Automatic adjustment of pacing amplitudes for improved patient safety and extended device longevity.

■ **DF4 connectors**

For simplified and shortened implantation procedures and less material in the device pocket.

■ **BIOTRONIK Home Monitoring®**

Enables wireless patient monitoring, including Home Monitoring-supported follow-ups with remote scheduling, for early detection of clinical and device-related events.

■ **Extended longevity**

Longer lasting device longevity through the use of energy-efficient technologies to reduce the number of device changes.

Ordering information

Model	Connectors	Volume/weight	Dimensions	Order number
Ilesto 7 VR-T	DF-1 (2x) IS-1 (1x)	33 cm ³ /80 g	65 mm x 55 mm x 11 mm	390082
Ilesto 7 VR-T	DF4 (1x)	31 cm ³ /80 g	65 mm x 52 mm x 11 mm	390088

CE₀₁₂₃

Cardiac Rhythm Management

Tachyarrhythmia Therapy

Single-chamber ICD

Ilesto 5 VR-T DX

MR Conditional single-chamber ICD with complete atrial diagnostics

ProMRI®



Product Highlights

■ **Complete atrial diagnostics**

Offers together with the Linx^{smart} ProMRI S DX lead the complete atrial information.

■ **ProMRI®**

Allows patients to undergo MR scanning under specific conditions.

■ **Small size**

Increases the patients' comfort through a reduced device thickness.

■ **SafeSync RF telemetry**

RF telemetry for wandless, time-saving and reliable data transmission at implantation and follow-up.

■ **Ventricular Capture Control**

Automatic adjustment of pacing amplitudes for improved patient safety and extended device longevity.

■ **Heart Failure Monitor**

Enables early detection of changes in patients' heart failure conditions by the continuous monitoring of clinical parameters.

■ **BIOTRONIK Home Monitoring®**

Enables wireless patient monitoring, including Home Monitoring-supported follow-ups with remote scheduling, for early detection of clinical and device-related events.

Ordering information

Model	Connectors	Volume/weight	Dimensions	Order number
Ilesto 5 VR-T DX	DF-1 [2x] IS-1 [2x]	33 cm ³ /81 g	65 mm x 55 mm x 11 mm	390122

CE₀₁₂₃

Cardiac Rhythm Management

Tachyarrhythmia Therapy

Single-chamber ICD

Ilesto 7 VR-T DX

MR Conditional single-chamber ICD with complete atrial diagnostics

ProMRI®



Product Highlights

■ **Complete atrial diagnostics**

Offers together with the Linx^{smart} ProMRI S DX lead the complete atrial information.

■ **ProMRI®**

Allows patients to undergo MR scanning under specific conditions.

■ **Small size**

Increases the patients' comfort through a reduced device thickness.

■ **Ventricular Capture Control**

Automatic adjustment of pacing amplitudes for improved patient safety and extended device longevity.

■ **Heart Failure Monitor**

Enables early detection of changes in patients' heart failure conditions by the continuous monitoring of clinical parameters.

■ **BIOTRONIK Home Monitoring®**

Enables wireless patient monitoring, including Home Monitoring-supported follow-ups with remote scheduling, for early detection of clinical and device-related events.

■ **Extended longevity**

Longer lasting device longevities through the use of energy-efficient technologies to reduce the number of device changes.

Ordering information

Model	Connectors	Volume/weight	Dimensions	Order number
Ilesto 7 VR-T DX	DF-1 (2x) IS-1 (2x)	33 cm ³ /81 g	65 mm x 55 mm x 11 mm	390094