



Cardiac Rhythm Management

Cardiac Resynchronization Therapy

CRT-D

## Ilesto 7 HF-T

MR Conditional three-chamber ICD

ProMRI®



### Product Highlights

■ **ProMRI®**

Allows patients to undergo MR scanning under specific conditions.

■ **Small size**

Increases the patients' comfort through a reduced device thickness.

■ **Capture Control**

Automatic adjustment of pacing amplitudes for improved patient safety and extended device longevity (in RA capture monitoring only).

■ **Heart Failure Monitor**

Enables early detection of changes in patients' heart failure conditions by the continuous monitoring of clinical parameters.

■ **DF4 connectors**

For simplified and shortened implantation procedures and less material in the device pocket.

■ **BIOTRONIK Home Monitoring®**

Enables wireless patient monitoring, including Home Monitoring-supported follow-ups with remote scheduling, for early detection of clinical and device-related events.

■ **Extended longevity**

Longer lasting device longevities through the use of energy-efficient technologies to reduce the number of device changes.

### Ordering information

Model	Connectors	Volume/weight	Dimensions	Order number
Ilesto 7 HF-T	DF-1 (2x) IS-1 (3x)	34 cm <sup>3</sup> /83 g	65 mm x 58.5 mm x 11 mm	390055
Ilesto 7 HF-T	DF4 (1x) IS-1 (2x)	33 cm <sup>3</sup> /82 g	65 mm x 56 mm x 11 mm	390061



Cardiac Rhythm Management

Cardiac Resynchronization Therapy

CRT-D

## Ilesto 5 HF-T

MR Conditional three-chamber ICD

ProMRI®



### Product Highlights

■ **ProMRI®**

Allows patients to undergo MR scanning under specific conditions.

■ **Small size**

Increases the patients' comfort through a reduced device thickness.

■ **SafeSync RF telemetry**

RF telemetry for wandless, time-saving and reliable data transmission at implantation and follow-up.

■ **Capture Control**

Automatic adjustment of pacing amplitudes for improved patient safety and extended device longevity (in RA capture monitoring only).

■ **Heart Failure Monitor**

Enables early detection of changes in patients' heart failure conditions by the continuous monitoring of clinical parameters.

■ **DF4 connectors**

For simplified and shortened implantation procedures and less material in the device pocket.

■ **BIOTRONIK Home Monitoring®**

Enables wireless patient monitoring, including Home Monitoring-supported follow-ups with remote scheduling, for early detection of clinical and device-related events.

### Ordering information

Model	Connectors	Volume/weight	Dimensions	Order number
Ilesto 5 HF-T	DF-1 (2x) IS-1 (3x)	34 cm <sup>3</sup> /83 g	65 mm x 58.5 mm x 11 mm	390110
Ilesto 5 HF-T	DF4 (1x) IS-1 (2x)	33 cm <sup>3</sup> /82 g	65 mm x 56 mm x 11 mm	390112